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RadiCS®
RadiCS® LE

Quality Control Software

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Önemli

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The screenshot shows the RadiCS LE software interface. On the left, the 'Device List' pane shows a tree structure under 'Computer' with 'Intel(R) HD Graphics 4600' and 'EIZO RX360' selected. The 'EIZO RX360' is expanded, showing options for 'DICOM', 'CAL1', 'CAL2', 'Custom', 'sRGB', and 'Text'. The 'DICOM' option is checked. On the right, the 'History List' pane shows a table of settings for the selected device.

Item	Value
CAL Switch Mode	DICOM
Calibration Target	Exp(2.2) (0.60cd/m ² -500.00cd/m ² Custom(x=0.2985, y=0.3104))
Hybrid Gamma PXL	☐ Enabled
Use/Comment	(undefined)
Backlight Meter	Insufficient amount of data
Backlight Status	Backlight is stable ✓

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Search condition

Monitor: ☒ Show only connected monitors

Keyword: AND OR

Result: ☐ Failed ☐ Passed ☐ Canceled ☐ Error ☐ Details / No Judgement / -

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Consistency Test

Luminance and Grayscale Check Failed.

100%

Detail

Luminance Measurement Result

L'max	174.72 cd/m ²	
L'min	0.26 cd/m ²	
Lamb	0.00 cd/m ²	

Result	Condition	Result
✓ Passed	L'max / L'min > 250	672.00
✓ Passed	L'max > 170.00 cd/m ²	174.72 cd/m ²
✓ Passed	Lamb < L'min / 1.5	0.17 cd/m ²
❌ Failed	ΔL'max < 10 % Base Value: 500.00 cd/m ²	-65.06 %

Grayscale Measurement Result

	Grayscale	Target Value	Measurement Value	Error Rate (%)
❌	0	0.60	0.26	-15.91
❌	15	1.54	0.60	-20.88
✓	30	3.10	1.07	-8.83
✓	45	5.46	1.80	-6.14
✓	60	8.85	2.83	0.78

Result	Condition	Max Error Rate
❌ Failed	Target Error Rate < 10% GSDF	37.39 %

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Assess whether the illuminance of EIZO RX360 is appropriate.

Room Category

You have selected Room Category RK1(<= 50lx). Please check if the current ambient illuminance is appropriate.

☐ Use the measurement value

Measurement Device

Serial Number (S/N)

Measurement Value lx

Click "Measure" to measure the illuminance with the monitor's illuminance sensor and automatically input the measurement value.
Input the value when using a measurement value from an external illuminance sensor.

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Assess whether the illuminance of EIZO RX360  DICOM is appropriate.

Application Category Application Category A(<= 50lx)

You have selected Application Category Application Category A(<= 50lx). Please check if the current illuminance is appropriate.

☐ Use the measurement value

Measurement Device

Serial Number (S/N)

Measurement Value lx Measure

Click "Measure" to measure the illuminance with the monitor's illuminance sensor and automatically input the measurement value.
Input the value when using a measurement value from an external illuminance sensor.

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Click "Proceed" to carry out the Calibration for EIZO RX360  DICOM.

Lamb ☐

During the Calibration, the following Lamb value is used. Please change or measure the value as required.

Measurement Device	LX-Can
Serial Number(S/N)	
Measurement Value	0.00 cd/m ²

[Measure](#)

* Turn off the monitor to measure the Lamb value manually.

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To retry the operation, select the monitor, and then click "Retry".

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Manually adjust the monitor brightness to within the Lmax target range by using the brightness settings.
Click "Start measurement" to measure monitor brightness.

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Home Device List **History List** Action Options

Search condition

Monitor ☒ Show only connected monitors Keyword AND OR

Result ☐ Failed ☐ Passed ☐ Canceled ☐ Error ☐ Details / No Judgement / -

Search results 13 Number of displays per page 100

Date	Job	Result	QC Guideline	Tester	Monitor	CAL Switch Mode
04/18/2019 13:21	Acceptance Test	Failed	DIN V 6868-57 Applicat...	RadiCS	EIZO RX360	DICOM
04/18/2019 13:21	Ambient luminance setting	Details	-	RadiCS	EIZO RX360	DICOM
04/18/2019 12:31	QC Guideline setting	Details	-	RadiCS	EIZO RX360	DICOM
04/18/2019 12:14	Baseline Value setting	Details	-	RadiCS	EIZO RX360	DICOM
04/18/2019 12:10	Calibration	Passed	-	RadiCS	EIZO RX360	DICOM
04/18/2019 11:56	Visual Check	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 11:56	Visual Check	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 11:46	Visual Check	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 11:44	Consistency Test(Biannual/Annual)	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 11:44	Consistency Test(Biannual/Annual)	Failed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 11:23	Acceptance Test	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM

History Import Bulk Test Report Generation

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The screenshot shows the RadiCS software interface. At the top, there's a navigation bar with 'Home', 'Device List', and 'History List' (selected). Below the navigation bar, there's a search condition section with a 'Monitor' dropdown, a checkbox for 'Show only connected monitors', a 'Keyword' input field, and radio buttons for 'AND' and 'OR'. Below this, there's a table of search results. The table has columns: Date, Job, Result, QC Guideline, Tester, Monitor, and CAL Switch Mode. The 'Delete' button is visible in the context menu for the selected test result.

Date	Job	Result	QC Guideline	Tester	Monitor	CAL Switch Mode
04/22/2019 14:39	Visual Check	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 13:21	Acceptance Test	Failed	DIN V 6868-57 Applicat...	RadiCS	EIZO RX360	DICOM
04/18/2019 13:21	Ambient luminance setting	Details	-	RadiCS	EIZO RX360	DICOM
04/18/2019 12:31	QC Guideline setting	Details	-	RadiCS	EIZO RX360	DICOM
04/18/2019 12:14	Baseline Value setting	Details	-	RadiCS	EIZO RX360	DICOM
04/18/2019 12:10	Calibration	Passed	-	RadiCS	EIZO RX360	DICOM
04/18/2019 11:56	Visual Check	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 11:56	Visual Check	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 11:46	Visual Check	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 11:44	Consistency Test(Biannual/Annual)	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 11:44	Consistency Test(Biannual/Annual)	Failed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM

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The screenshot shows the RadiCS software interface. At the top, there are tabs for 'Home', 'Device List', and 'History List'. Below the tabs is a search bar with a 'Search condition' section. The 'Monitor' dropdown is set to 'All'. The 'Show only connected monitors' checkbox is checked. The 'Keyword' field is empty. The 'Result' dropdown is set to 'Failed'. The 'Search results' section shows 14 results. The table below lists the test results:

Date	Job	Result	QC Guideline	Tester	Monitor	CAL Switch Mode
04/22/2019 14:39	Visual Check	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 13:21	Acceptance Test	Failed	DIN V 6868-57 Applicat...	RadiCS	EIZO RX360	DICOM
04/18/2019 13:21	Ambient luminance setting	Details	-	RadiCS	EIZO RX360	DICOM
04/18/2019 12:31	QC Guideline setting	Details	-	RadiCS	EIZO RX360	DICOM
04/18/2019 12:14	Baseline Value setting	Details	-	RadiCS	EIZO RX360	DICOM
04/18/2019 12:10	Calibration	Passed	-	RadiCS	EIZO RX360	DICOM
04/18/2019 11:56	Visual Check	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 11:56	Visual Check	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 11:46	Visual Check	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 11:44	Consistency Test(Biannual/Annual)	Passed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM
04/18/2019 11:44	Consistency Test(Biannual/Annual)	Failed	JESRA Grade 1A	RadiCS	EIZO RX360	DICOM

At the bottom of the table, there are buttons for 'History Import' and 'Bulk Test Report Generation'. A 'Show report' button is also visible over the table.

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

The dialog box is titled "Illuminance Sensor Correlation" and contains the following text: "Execute the Illuminance Sensor Correlation for EIZO RX350-AR [monitor icon]. Enter the measured value using an external illuminance sensor, then click 'Proceed'." Below the text is a checkbox labeled "Reset illuminance sensor correlation". Under the checkbox is a text input field with the unit "lx" to its right. At the bottom right are two buttons: "Cancel" and "Proceed".

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The window shows the "Illuminance Sensor Correlation Result" with a progress bar at the top indicating three steps: "1 Execution options", "2 Proceed Illuminance Sensor Correlation" (which is the active step), and "3 Finish". Below the progress bar is a table with the following data:


Monitor	Remarks
EIZO RX350-AR [monitor icon]	Finished.

At the bottom left is a "Cancel" button and at the bottom right is a "Finish" button.

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RadiCS						
Execution timing	Monitor	CAL Swit...	Job	Tester	Duration	Status
08/01/2022 00:00	EIZO MX216	DICOM	Consistency Test	RadiCS(Scheduled)	-	Unexecuted

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
The screenshot shows the RadiCS software interface. The top menu bar includes 'Home', 'Device List', 'History List', 'Action', and 'Options'. The 'Options' menu is open, showing a sub-menu with 'Enable Hide-and-Seek' checked. Below this, there is a table with columns: Monitor, Display Position, Switch position, and Hotkey. The table contains one row with the following data: Monitor: EIZO RX360, Display Position: Upper Right, Switch position: Top center edge + Top right edge + ..., Hotkey: None. Below the table are buttons for 'Add', 'Delete', and 'Edit'. At the bottom right of the 'Options' menu are buttons for 'Undo' and 'Save'.

Monitor	Display Position	Switch position	Hotkey
EIZO RX360	Upper Right	Top center edge + Top right edge + ...	None

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1. Select the monitor to enable the PinP sub-window

EIZO RX360

2. Configure the PinP sub-window position

Upper Right

Offset

Horizontal (H) 0 px

Vertical (V) 0 px

3. Configure the mouse/hotkey operation of the PinP sub-window

☒ Mouse operation

Detection position Click the dotted line shown in the figure to configure the detection position.

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
RadiCS Version 5.1.1 About RadiCS

Home Device List History List Action Options

Monitor input signals can be switched according to the hotkey or Switch-and-Go action.

☒ Enable Signal Switch

1. Select monitors and sets of input signals.

Select an input signal available on the monitor.

Monitor	Input Signal 1	Input Signal 2
☒ EIZO RX360	DP1	DSUB1
☐ EIZO RX360	DP1	DSUB1

2. Select a timing to switch the Input Signal.

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Please enter the hotkey directly.
Functions with * cannot use a duplicate hotkey.

Function	Monitor	Hotkey
Hide-and-Seek		None
Switch-and-Go		None
Point-and-Focus *		None
Manual Mode Switch *		None
Signal Switch		None
Move to home position *		None
Painter Position Indication *		None

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Item	Value
Manufacturer	Intel Corporation
Serial Number	(undefined)
Driver	igfx
Driver Version	9.18.10.3204
Installed on	09/05/2018

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The screenshot shows the RadiCS software interface. On the left, the 'Device List' tab is active, showing a tree structure under 'Computer' with 'Intel(R) HD Graphics 4600' and two EIZO monitors: 'EIZO RX360' and 'EIZO RX1270'. The 'EIZO RX1270' monitor is selected. On the right, a detailed view of the selected monitor is shown, including a table of properties.

Item	Value
Asset Number	[undefined]
Usage Time (Daily Average)	8H (-)
Installed on	10/17/2019
Connection	USB
Luminance Sensor	Integrated Front Sensor
Presence Sensor	-
Illuminance Sensor	Yes
Key Lock	OFF
Size in inches	30.9
Resolution	4200x2800 @ 29Hz
Monitor Type	Color (Hardware Calibration)
UDI	-
RadiLight Area	RadiLight Area: ON, Brightness: 5

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The screenshot shows the RadiCS software interface. The 'Options' menu is open, and the 'Change...' button is visible next to the Password field. The interface includes a sidebar with various settings categories like General, Registration Information, Schedule, Sensor, User Mode, History, Ambient Light Watchdog, and MAC Address Clone. The main area displays settings for Password, Illuminance, SelfQC History, Tester, Monitor Detection, Language, Loglevel, and Remote Setting.

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The screenshot shows the 'Parola Ayarla' (Set Password) dialog box. It contains three input fields: 'Current Password', 'New Password', and 'Type New Password'. There are 'Cancel' and 'OK' buttons at the bottom right.

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9 Information

This chapter provides the following information:

- Notes concerning the monitor quality control standards (QC guidelines) used by RadiCS.
- Precautions for setting up a test in RadiCS according to each monitor quality control standards (QC guidelines).

9.1 Description of Standards

9.1.1 Quality Control Standards for Digital Imaging for Medical Display Monitors (Monitor Quality Control Standards)

IEC 62563-2: 2021

"Medical electrical equipment - Medical image display systems - Part 2: Acceptance and constancy tests for medical image displays" issued by the International Electrotechnical Commission. This standard uses the evaluation method of IEC 62563-1 to specify test criteria, frequency, category classification, etc.

Not	
• "IEC 62563-2" in RadiCS includes the following.	
Standard / Guideline references	QC guideline (Abbreviation)
Category I-A	IEC 62563-2 Category I-A
Category I-B	IEC 62563-2 Category I-B
Category II ^{*1}	IEC 62563-2 Category II for Diagnosis
	IEC 62563-2 Category II for Viewing
^{*1} Category II is divided into two categories in RadiCS because the evaluation contents and judgment criteria are different for diagnostic and viewing purposes.	

AAPM On-line Report No. 03: 2005

"Assessment of Display Performance for Medical Imaging Systems" formulated by Task Group (TG) 18 of American Association of Physicists in Medicine. It defines consistency tests and acceptance tests for monitors. Monitors are classified into "Primary" and "Secondary" depending on the intended use.

Not
• "AAPM" used in RadiCS means "AAPM On-line Report No. 03".

ACR-AAPM-SIIM Practice Guideline for Determinants of Image Quality in Digital Mammography: 2012

This guideline was formulated collaboratively by specialists in mammography and medical physics who represent the American College of Radiology (ACR), American Association of Physicists in Medicine (AAPM), and Society for Imaging Informatics in Medicine (SIIM). The Mammography Quality Standards Act (MQSA) obliges the quality control for mammography diagnostic equipments in the United States. This Act, which went into effect in 1992, is aimed at film based analog systems, and is being revised for digital systems that become popular recently. This guideline is positioned as one of proposals by ACR for such rework. The section on monitors covers diagnostic (Primary) use. It does not cover the concepts of acceptance tests or consistency tests. This was revised in 2012.

Not

- RadiCS with "ACR" indicates that it has been tested with additional quality control elements based on the ACR-AAPM-SIIM Practice Guideline for Determinants of Image Quality in Digital Mammography (hereinafter referred to as ACR Mammo) (the evaluation item and standard are selected from the ACR-AAPM-SIIM Practice Guideline for Determinants of Image Quality in Digital Mammography: 2012 (hereinafter referred to as the Technical Standard) and AAPM On-line Report No. 03:2005).

New York State Department of Health Bureau of Environmental Radiation Protection Guide for Radiation Safety / Quality Assurance Program Primary Diagnostic Monitors

The guidelines describe the types and extension of information and criteria used by the New York State Department of Health Bureau of Environmental Radiation Protection to evaluate Primary Diagnostic Monitor (PDM) in facilities as a part of the radiation safety and quality assurance program.

Not

- Term "NYS PDM-****" in RadiCS refers to "New York State Department of Health Bureau of Environmental Radiation Protection Guide for Radiation Safety/Quality Assurance Program Primary Diagnostic Monitors". In RadiCS, contents are added by referring partially to AAPM On-line Report No. 03.

Standard / Guideline references	QC guideline (Abbreviation)
Not for mammography	NYS PDM – Diagnostic
For mammography	NYC PDM – Clinical sites

NYC Quality Assurance Guidelines for Primary Diagnostic Monitors: 2015

Refers to the "Guidance related to quality assurance for Primary Diagnostic Monitor (PDM)" based on the health regulations of New York city provided by the New York City Health Department's Office of Radiological Health.

Not	
The term "NYC PDM-****" in RadiCS refers to "NYC Quality Assurance Guidelines for Primary Diagnostic Monitors: 2015". In RadiCS, contents are added by referring partially to AAPM On-line Report No. 03.	
Standard / Guideline references	QC guideline (Abbreviation)
For hospitals, medical centers, imaging centers, radiologist offices	NYC PDM – Hospitals
For all other clinical sites, including chiropractic offices, medical doctor offices, orthopedic offices	NYC PDM – Clinical sites
For mammography facilities	NYC PDM – Mammography

ONR 195240-20: 2017

"Image Quality Assurance in X-ray Diagnosis - Part 20: Acceptance test and consistency test for image display devices" formulated by the Austrian Standards Institute. This standard is based on German DIN 6868-157 and QS-RL standards, with the Institute's own judgment and interpretation added to the compilation. Compared with the 2008 edition, parts of test patterns, evaluation methods, judgment standards, etc. to be used have been modified in the new edition.

Not	
The term "ONR 195240-20 ***" in RadiCS refers to "Image Quality Assurance in X-ray Diagnosis - Part 20: Acceptance test and consistency test for image display devices: 2017".	
Standard / Guideline references	QC guideline (Abbreviation)
Mammography: Application Category A	ONR 195240-20 Application Category A Mammo
Application Category A	ONR 195240-20 Application Category A
In dentistry: Application Category B	ONR 195240-20 Application Category B Dentistry
Application Category B	ONR 195240-20 Application Category B

DIN 6868-157: 2022

"Image quality assurance in diagnostic X-ray – Part 157: X-ray Ordinance Acceptance and Consistency Tests of image display systems in their environment" formulated by the German Institute for Standardization (Deutsches Institut für Normung e.V). The standard is intended to replace the preceding DIN V 6868-57 standard that defines acceptance testing and the corresponding chapters of QS-RL and PAS1054 (see below) that specifies criteria by body part and capture method, consistency test items, and frequencies. Conformance to the international standard is also one of the reasons of revision and many of the evaluation methods and test patterns specified in IEC 62563-1 (or DIN EN 62563-1) have been adapted. There are also original approaches such as definition of room category and setting down of upper limit of illuminance according to the application. RadiCS reflects relevant items according to "QS-RL Rundschreiben (TOP C 04 der 74. Sitzung des LA RöV im Mai 2015, TOP C 07 der 75. Sitzung des LA RöV im November 2015)".

Not

- "DIN 6868-157" shown in RadiCS includes the followings.

Standard / Guideline references	QC guideline (Abbreviation)
DIN 6868-157 I. Mammography	DIN 6868-157 I. Mammography
DIN 6868-157 II. Mammographic stereotaxy	DIN 6868-157 II. Mammographic stereotaxy
DIN 6868-157 III. Projection radiography (thorax, skeleton, abdomen)	DIN 6868-157 III. Projection radiography
DIN 6868-157 IV. Fluoroscopy, all applications	DIN 6868-157 IV. Fluoroscopy, all applications
DIN 6868-157 V. Computed tomography	DIN 6868-157 V. Computed tomography
DIN 6868-157 VI. Digital volume tomography(dental), intraoral X-ray diagnostics with dental tubehead, panoramic radiographs, cephalometric radiographs of the skull, Dental radiographs of a skull overview, Hand radiographs for skeletal growth determination	DIN 6868-157 VI. Digital volume tomography (dental) etc. in RK 5
DIN 6868-157 VII. Intraoral X-ray diagnostics with dental tubehead, panoramic radiographs, cephalometric radiographs of the skull, Dental radiographs of a skull overview, Hand radiographs for skeletal growth determination (The interval of the measuring tests can be extended to five years on the condition that the requirements specified in TOP C 07 der 75. Sitzung des LA RöV are satisfied.)	DIN 6868-157 VI. Dental X-ray equipment etc. in RK 5 (five-year interval)
DIN 6868-157 VII. Intraoral X-ray diagnostics with dental tubehead, panoramic radiographs, cephalometric radiographs of the skull, Dental radiographs of a skull overview, Hand radiographs for skeletal growth determination	DIN 6868-157 VII. Intraoral X-ray diagnostics (dental) etc. in RK6
DIN 6868-157 VIII. Viewing	DIN 6868-157 VIII. Viewing

DIN V 6868-57: 2001

"Image Quality Assurance in X-ray Diagnosis - Part 57: Acceptance test for image display devices" formulated by the German Institute for Standardization (Deutsches Institut für Normung e.V.). Image display devices are divided into three categories. "Application Category A" includes image display devices used for the diagnosis of images of high spatial and contrast resolution. "Application Category B" includes image display devices for diagnosis which are not classified in "Application Category A" and image display devices for image viewing.

Quality Control Manual for Digital Mammography: 2017

A quality control manual for digital mammography systems written by the Japan Central Institute on Quality Assurance of Breast Cancer Screening, a nonprofit organization, in Japan. This NPO studies and manages quality control of mammography.

Not

- "DMG QC Manual" or "DMG QCM" in RadiCS refers to "Quality Control Manual for Digital Mammography". Note that "Regular Control Point" or "Daily Control Point" written in the DMG QCM is expressed as "Consistency Test" or "Visual Check" on RadiCS.

European Guidelines for Quality Assurance in Breast Cancer Screening and Diagnosis Fourth Edition - Supplements: 2013

This guideline was issued by the European Commission in cooperation with EUREF (European Reference Organization for Quality Assured Breast Screening and Diagnostic Services), EBCN (European Breast Cancer Network), and EUSOMA (European Society of Mastology). It applies to mammography systems as a whole and chapter 2 deals with monitors. Supplements were added in 2013. Different conditions are set for monitors for diagnostic and for reference use.

Not

- "EUREF" written on RadiCS means "European Guidelines for Quality Assurance in Breast Cancer Screening and Diagnosis Fourth Edition - Supplements".

JESRA TR-0049⁻²⁰²⁴

It refers to the "Guidelines for Acceptance and Consistency Tests of Medical Imaging Display Systems (JIS T 62563-2)" prepared by Japan Medical Imaging and Radiological Systems Industries Association (JIRA). It was established in 2024 to align with JIS and replace the previously widely used JESRA X-0093. This standard uses the evaluation method of JIS T 62563-1 to specify the test criteria, frequency, category classifications, and other aspects defined in JIS T 62563-2. Compliance with this standard also signifies compliance with JIS T 62563-2.

In the appendix to JESRA TR-0049, a category called "Category II (for diagnosis plus)," which emphasizes consistency between Category II (for diagnosis) and management grade 1B of JESRA X-0093, and visual tests for each date of use are posted as reference information. RadiCS also implements these features.

Not

- In RadiCS, any reference to "JESRA TR-0049" refers to the following.

Standard / Guideline references	QC guideline (Abbreviation)
Category I-A	JESRA TR-0049 Category I-A
Category I-B	JESRA TR-0049 Category I-B
Category II (for diagnosis plus)	JESRA TR-0049 Category II (for diagnosis plus)
Category II (for diagnosis)	JESRA TR-0049 Category II (for diagnosis)
Category II (for reference)	JESRA TR-0049 Category II (for reference)

JESRA X-0093*B⁻²⁰¹⁷: 2017

"Quality Assurance (QA) Guideline for Medical Imaging Display Systems" prepared by Japan Medical Imaging and Radiological Systems Industries Association (JIRA). It was published in 2005 and revised in 2010 and 2017. This guideline specifies the acceptance tests and consistency tests. Also, in this guideline, the organization can omit the acceptance test by substituting it with the shipment test reports provided by manufacturers. In the 2017 revision, the previous "Grade 1" was changed to "Grade 1B", and the new "Grade 1A" was added as the higher-level judgment criteria. The organization must judge which grade level is to be used for management depending on the intended use.

IPEM Report 91: 2005

"Recommended Standards for the Routine Performance Testing of Diagnostic X-ray Imaging Systems" formulated by Institute of Physics and Engineering in Medicine in the UK. It applies to diagnostic X-ray imaging systems as a whole including image display devices but does not include MR or ultrasonic systems. The items related to monitors were added when this standard was revised from Report 77. It mainly defines consistency tests.

Not

- "IPEM" used in RadiCS means "IPEM Report 91".

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"Guideline for implementing quality assurance of the X-ray systems for diagnostic and medical treatment purposes according to chapters 16 and 17 of the X-ray Ordinance". This defines the details of the quality assurance of general X-ray systems obliged by the X-ray Ordinance (for diagnostics: chapter 16, for medical treatment: chapter 17). DIN V 6868-57 is supposed to be referred on basic test methods for diagnostic image display devices. Limiting values such as the minimum value of the maximum luminance and the items/frequency of the consistency test are added to the contents of DIN V6868-57 that defines only the acceptance test. Although the classification of image display devices conforms to DIN V 6868-57 (Category A, B), stricter criteria are established for mammography equipments by reference to PAS1054 "Requirements and testing of digital mammographic X-ray equipment", which is the standard issued by the German Institute for Standardization.

Not

- "QS-RL" used in RadiCS means "Qualitätssicherungs-Richtlinie: 2007". "Application Category A Mammo" means PAS1054 is also complied with.

9.1.2 Other Standards

DICOM PS 3.14: 2000

"Digital Imaging and Communications in Medicine (DICOM) Part 14: Grayscale Standard Display Function" formulated by NEMA (National Electrical Manufacturers Association) in the US. It defines the grayscale characteristics to be equipped in films and monitors for the display of grayscale images as GSDF: Grayscale Standard Display Function. More details on the evaluation of compliance for this standard are specified in other policies and standards, such as AAPM On-line Report No. 03.

Not

- "DICOM Part 14 GSDF" used in RadiCS means "The grayscale standard display function defined in DICOM PS 3.14".

CIE Pub.15.2: 1986

"Colorimetry, Second Edition" published by Commission Internationale de l' Eclairage. It recommends CIELAB($L^*a^*b^*$) and CIELUV($L^*u^*v^*$) that are uniform color spaces and uses color difference formulas to evaluate the difference of two colors quantitatively.

Not

- "CIE" used in RadiCS means "Display formulas with L^* formula".

SMPTE RP133: 1991

"Specifications for Medical Diagnostic Imaging Test Pattern for Television Monitors and Hard-Copy Recording Cameras" proposed by Society of Motion Picture and Television Engineers in the US.

Not

- "SMPTE" used in RadiCS means "Test patterns created in reference to SMPTE RP133 specifications".

Basic QC, Basic Mammo QC, Basic Mammo QC for Remote, Basic QC Primary, Basic QC Primary for Remote, Basic QC Secondary, Basic QC Secondary for Remote, Pathology350, Pathology450

The setting specific to RadiCS used for monitor management that does not comply with standards or guidelines established in each country.

9.2 RadiCS Software

9.2.1 Prerequisite

RadiCS software

We have long developed monitors. With those skills, knowledge and measuring data, we have developed RadiCS for users of digital imaging for medical diagnosis to manage the quality of monitors efficiently according to our interpretation of the quality control standard for each digital imaging for medical monitor.

Each digital imaging for medical monitor evaluation standard defines the change of clinical image use and monitor luminance, as well as measurement devices. Having only RadiCS will not meet all the conditions. Read thorough the related standards and test each item according to the conditions.

A setting value for each standard can be changed and testing conditions can be set with several standards.

To maintain and manage image quality according to the standards and the situation, follow the monitor quality control standards and use RadiCS.

Monitor judgment by RadiCS is not to ensure each monitor quality control standard.

This product includes open source software.

If the open source software contains a product for which usage is granted under a GPL (GNU GENERAL PUBLIC LICENSE) license, EIZO Corporation will, in line with the GPL usage license conditions, provide the source code for corresponding GPL software via a medium, such as CD-ROM, at a cost to individuals and organizations who make contact via the following contact information for a minimum period of three years after purchase of the product.

We will also provide the source code for corresponding LGPL (GNU LESSER GENERAL PUBLIC LICENSE) software of products that include LGPL software licensed under the LGPL in the same manner as stated above.

Contact information

www.eizoglobal.com/contact/index.html

*Contact your local EIZO representative.

Except for open source software licensed under GPL, LGPL or other licenses, any transferring, copying, reverse assembly, reverse compiling or reverse engineering of any software included with this product is strictly prohibited. Further, exporting of any software included with this product in violation of applicable export laws is strictly prohibited.

9.2.2 Correlation Between RadiCS and Monitor Quality Control Standards

The RadiCS software interprets and supports each monitor quality control standard as described below. Use this information when setting up tests in RadiCS.

IEC 62563-2

RadiCS Setup

	Acceptance Test			
	Category I-A	Category I-B	Category II for Diagnosis	Category II for Viewing
Pattern Check (Used pattern)	TG18-OIQ TG18-MP TG18-UN80 TG18-UN10	TG18-OIQ TG18-MP TG18-UN80 TG18-UN10	TG18-OIQ TG18-MP TG18-UN80	TG18-OIQ TG18-MP TG18-UN80
Luminance Check	$L'_{max} > 450 \text{ cd/m}^2$ $L'_{max} / L'_{min} > 350$ $L_{amb} < L_{min} / 0.67$	$L'_{max} > 350 \text{ cd/m}^2$ $L'_{max} / L'_{min} > 250$ $L_{amb} < L_{min} / 0.67$	$L'_{max} > 150 \text{ cd/m}^2$ $L'_{max} / L'_{min} > 100$	$L'_{max} > 150 \text{ cd/m}^2$ $L'_{max} / L'_{min} > 100$
Grayscale Check	Target error rate < 10 % of GSDF Grayscale chromaticity delta $u'v' < 0.010$ (5.00 cd/m^2 or more)	Target error rate < 10 % of GSDF Grayscale chromaticity delta $u'v' < 0.010$ (5.00 cd/m^2 or more)	Target error rate < 20 % of GSDF Grayscale chromaticity delta $u'v' < 0.015$ (5.00 cd/m^2 or more)	Target error rate < 20 % of GSDF
Uniformity Check	Grayscale 204 $(L_{max}-L_{min}) / (L_{max}+L_{min}) \times 200 < 20 \%$ Grayscale 204 $\Delta u'v' < 0.010$	Grayscale 204 $(L_{max}-L_{min}) / (L_{max}+L_{min}) \times 200 < 20 \%$ Grayscale 204 $\Delta u'v' < 0.010$	Grayscale 204 $(L_{max}-L_{min}) / (L_{max}+L_{min}) \times 200 < 30 \%$ Grayscale 204 $\Delta u'v' < 0.015$	Grayscale 204 $(L_{max}-L_{min}) / (L_{max}+L_{min}) \times 200 < 30 \%$
Multi-monitor	$\Delta L'_{max} < 10 \%$ Grayscale 204 $\Delta u'v' < 0.010$	$\Delta L'_{max} < 10 \%$ Grayscale 204 $\Delta u'v' < 0.010$	$\Delta L'_{max} < 20 \%$ Grayscale 204 $\Delta u'v' < 0.015$	$\Delta L'_{max} < 20 \%$
Resolution	$\geq 2048 \times \geq 2048$	$\geq 1024 \times \geq 1024$	-	-

	Consistency Test			
	Category I-A	Category I-B	Category II for Diagnosis	Category II for Viewing
Pattern Check (Used pattern)	TG18-OIQ TG18-MP TG18-UN80 TG18-UN10	TG18-OIQ TG18-MP TG18-UN80 TG18-UN10	TG18-OIQ TG18-MP TG18-UN80	TG18-OIQ TG18-MP TG18-UN80
Luminance Check	$L'_{max} > 450 \text{ cd/m}^2$ $L'_{max} / L'_{min} > 350$ $L_{amb} / L_{min} < 0.67$	$L'_{max} > 350 \text{ cd/m}^2$ $L'_{max} / L'_{min} > 250$ $L_{amb} / L_{min} < 0.67$	$L'_{max} > 150 \text{ cd/m}^2$ $L'_{max} / L'_{min} > 100$	$L'_{max} > 150 \text{ cd/m}^2$ $L'_{max} / L'_{min} > 100$
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 10 % of GSDF	Target error rate < 20 % of GSDF	Target error rate < 20 % of GSDF
Uniformity Check	-	-	-	-
Multi-monitor	$\Delta L'_{max} < 10 \%$	$\Delta L'_{max} < 10 \%$	$\Delta L'_{max} < 20 \%$	$\Delta L'_{max} < 20 \%$
Resolution	$\geq 2048 \times \geq 2048$	$\geq 1024 \times \geq 1024$	-	-

IEC 62563-2: 2021 and RadiCS

Pattern Check

RadiCS prepares the patterns based on check results for respective compatible resolutions.

Luminance Check

The standard includes an equality sign in each judgment condition but RadiCS does not include an equality sign.

The " L_{amb}/L'_{min} (a) relationship <0.6 " equation has been changed to " $L_{amb} < L'_{min}/0.67$ " to determine the ambient luminance.

Grayscale Check

The standard includes an equality sign in each judgment condition but RadiCS does not include an equality sign.

In RadiCS, "target error rate < 10 or 20 % of GSDF" indicates a contrast response test, which measures 18 points. Measured values of less than 5.00 cd/m^2 are not used to determine "Grayscale chromaticity $\Delta u'v'$ ".

Uniformity Check

The standard includes an equality sign in each judgment condition but RadiCS does not include an equality sign.

It describes how to use the TG18-UNL80 pattern, but RadiCS displays a 10 % display area of the window at grayscale 204 in the middle and corner of the screen, and measures the center of the window.

Sensors

Noncontact and contact measurement devices can be used in IEC 62563-2.

Multi-monitor

The standard includes multi-monitor judgment and includes an equality sign, but RadiCS does not include an equality sign.

Cautions

Although Category II is not classified in the standard, RadiCS divides it into two categories for convenience, since the evaluation contents/judgment criteria differ between diagnostic and viewing use. Note that Category III in the standard is not implemented in RadiCS.

AAPM**RadiCS Setup**

	Acceptance Test	
	Primary	Secondary
Pattern Check (Used pattern)	Black TG18-QC TG18-AD TG18-UN80 TG18-AFC TG18-CT White	Black TG18-QC TG18-AD TG18-UN80 TG18-AFC TG18-CT White
Luminance Check	$L'_{\max} / L'_{\min} > 250$ $L'_{\max} > 170 \text{ cd/m}^2$ $\Delta L'_{\max} < 10 \%^{*1}$	$L'_{\max} / L'_{\min} > 100$ $L'_{\max} > 100 \text{ cd/m}^2$ $\Delta L'_{\max} < 10 \%^{*1}$
Grayscale Check	Target error rate $< 10 \%$ of GSDF	Target error rate $< 20 \%$ of GSDF
Uniformity Check	Grayscale: 204, 26 ^{*2} Grayscale: 204 $\Delta u'v' < 0.010$	Grayscale: 204, 26 ^{*2}
Multi-monitor	$\Delta L'_{\max} < 10 \%$ between multiple monitors Grayscale 204 Mean value between multiple monitors $\Delta u'v' < 0.010$	$\Delta L'_{\max} < 10 \%$ between multiple monitors

	Consistency Test	
	Primary	Secondary
Pattern Check (Used pattern)	TG18-QC TG18-AD TG18-UN80 TG18-AFC Black White	TG18-QC TG18-AD TG18-UN80 TG18-AFC Black White
Luminance Check	$L'_{\max} / L'_{\min} > 250$ $L'_{\max} > 170 \text{ cd/m}^2$ $\Delta L'_{\max} < 10 \%^{*1}$	$L'_{\max} / L'_{\min} > 100$ $L'_{\max} > 100 \text{ cd/m}^2$ $\Delta L'_{\max} < 10 \%^{*1}$
Grayscale Check	Target error rate $< 10 \%$ of GSDF	Target error rate $< 20 \%$ of GSDF
Uniformity Check	Grayscale: 204, 26 ^{*2}	Grayscale: 204, 26 ^{*2}
Multi-monitor	$\Delta L'_{\max} < 10 \%$ between multiple monitors	$\Delta L'_{\max} < 10 \%$ between multiple monitors

*1 $L_{\text{amb}} < L_{\min} / 1.5$ *2 $(L_{\max} - L_{\min}) / (L_{\max} + L_{\min}) \times 200 < 30 \%$

Correlation between AAPM and RadiCS

Pattern Check

A test pattern given in AAPM cannot be applied to a monitor whose screen aspect ratio is not 1:1 without modification, since AAPM (or the test pattern) uses an aspect ratio of 1:1. Therefore, RadiCS checks a monitor being tested, and determines and generates an appropriate test pattern for each resolution supported by the monitor.

TG18-QC	Equivalent to the pattern with the same name in the standard. The pattern is scaled in accordance with the screen resolution.
TG18-AD	
TG18-AFC	
TG18-CT	
TG18-UN80	Grayscale 204 white patterns. The same pattern of AAPM has a square frame but RadiCS does not have any because it does not need to be visible.

Luminance Check

AAPM except for $L_{amb} < L_{min}$ includes an equality sign in each judgment condition but RadiCS does not.

The calibration setup, L_{max} value will be input in the $\Delta L'_{max}$ baseline value as an initial setup when performing a tasksetup.

L'_{max}/L'_{min} means AAPM $LR' (= (L_{max} + L_{amb}) / (L_{min} + L_{amb}))$.

Grayscale Check

AAPM includes an equality sign but RadiCS doesn't because of the target error rate is $< 10\%$ of GSDF. This is a judgment condition for DICOM Part 14 GSDF.

The number of grayscale measuring points is fixed at 18 and is unchangeable.

The measurement result is 17 points because it is expressed as $(JND_{n+1} - JND_n) / 2$.

Uniformity Check

AAPM includes an equality sign in each judgment condition but RadiCS does not.

AAPM uses TG18-UN80 and TG18-UN10 patterns in measurement, but these patterns cannot be applied to a monitor whose screen aspect ratio is not 1:1 without modification, since they use an aspect ratio of 1:1. Instead, RadiCS displays grayscale 204 and grayscale 26 windows equivalent to 10 % of the display area in the center of the screen and in the corners, and measures the center portion of each window.

Sensors

Noncontact and contact measurement devices are available in AAPM.

Multi-monitor

AAPM includes a determination for multiple monitors, but by default RadiCS is set not to make such a determination. If necessary, make settings as indicated in the table above. AAPM includes an equality sign but RadiCS does not.

Cautions

AAPM consistency testing has three types: tests that monitor users perform daily, tests that medical physicists perform or QC (quality control) technologists perform under their instructions monthly / quarterly, and tests that medical physicists perform annually. RadiCS is mainly intended for consistency testing of the second type, but pattern checks can be performed for all three types of testing.

AAPM has an item to measure geometrical distortion but RadiForce series monitors do not need to be measured because it meets the requested specification.

However, non-RadiForce monitors may be used. Therefore, the pattern check has patterns and checkpoints for geometrical distortion.

ACR**RadiCS Setup**

	Acceptance Test	Consistency Test
Pattern Check (Used pattern)	Black TG18-QC TG18-AD TG18-UN80 TG18-AFC TG18-CT White	Black TG18-QC TG18-AD TG18-UN80 TG18-AFC White
Luminance Check	$L'max / L'min > 250$ $L'max > 420 \text{ cd/m}^2$ $L'min > 1.2 \text{ cd/m}^2$ $Lamb < Lmin / 4$	$L'max / L'min > 250$ $L'max > 420 \text{ cd/m}^2$ $L'min > 1.2 \text{ cd/m}^2$ $Lamb < Lmin / 4$
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 10 % of GSDF
Uniformity Check	Grayscale: 204, 26 $(Lmax-Lmin) / (Lmax+Lmin) \times 200 < 30 \%$ Grayscale: 204 $\Delta u'v' < 0.010$	Grayscale: 204, 26 $(Lmax-Lmin) / (Lmax+Lmin) \times 200 < 30 \%$ Grayscale: 204 $\Delta u'v' < 0.010$
Multi-monitor	Grayscale: 204 $\Delta u'v' < 0.010$	Grayscale: 204 $\Delta u'v' < 0.010$

Correlation between ACR and RadiCS**Pattern Check**

The test patterns are not introduced specifically in ACR Mammo. The same check method as AAPM is applied to RadiCS. See the AAPM item for details of the correlation with RadiCS.

Luminance Check

For ACR Mammo, only " $L'max \geq 400 \text{ cd/m}^2$ (recommendation: 450 cd/m^2)" is displayed. For the Technical Standard, " $L'max \geq 420 \text{ cd/m}^2$ " is specified for mammography, so 420 cd/m^2 is used. In addition, other judgment standards specified by the Technical Standard are also used. The judgment conditions include an equality sign but RadiCS does not.

Grayscale Check

GSDF is recommended for ACR Mammo, but there is no judgment standard. For reference values, the values for AAPM and the Technical Standard are used. These include an equality sign but RadiCS does not because the target error rate is < 10 % of GSDF. This is a judgment condition for DICOM Part 14 GSDF.

The number of grayscale measuring points is fixed at 18 and is unchangeable. The measurement result is 17 points because it is expressed as $(JND_{n+1} - JND_n)/2$.

Uniformity Check

For ACR Mammo, the uniformity of the luminance and chromaticity is not specified. The uniformity needs to be confirmed, so conditions for RadiCS include uniformity judgment for luminance and chromaticity. The content is the same as that for AAPM. For details on the correlation with RadiCS, see the AAPM section.

Sensors

ACR Mammo contains nothing in particular about sensors or measurement devices. Since this standard was compiled using AAPM as a reference, sensors are handled in the same manner as AAPM.

Multi-monitor

For ACR Mammo, there is no multi-monitor judgment. By default, RadiCS does not perform judgment. If necessary, make settings as indicated in the table above.

Cautions

ACR Mammo is an educational tool to supply physicians, technicians, and physicists with extensive knowledge related to digital mammography image quality. It is not an implementation standard, a list of essential requirements, or a quality control standard. For this reason it does not cover the concepts of acceptance tests or consistency tests. However, we, who have agreed to the ACR policy, suggest support for the deficiencies in quality control with reference to the AAPM and the Technical Standard stated in ACR Mammo to achieve more practical operation.

NYS-PDM**RadiCS Setup**

	Acceptance Test / Consistency Test [Annually]	
	NYS PDM – Diagnostic	NYS PDM – Mammography
Pattern Check (Used pattern)	-	-
Luminance Check	L'max / L'min > 170 L'max > 171 cd/m ² Lamb < Lmin / 1.5	L'max / L'min > 250 L'max > 250 cd/m ² Lamb < Lmin / 1.5
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 10 % of GSDF
Uniformity Check	Grayscale: 204, 26 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 30 %	Grayscale: 204, 26 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 30 %
Multi-monitor	-	-

	Consistency Test [Bi-Weekly]	
	NYS PDM – Diagnostic	NYS PDM – Mammography
Pattern Check (Used pattern)	Black SMPTE Shades of RGB White	Black SMPTE Shades of RGB White
Luminance Check	-	-
Grayscale Check	-	-
Uniformity Check	-	-
Multi-monitor	-	-

	Consistency Test [Quarterly]	
	NYS PDM – Diagnostic	NYS PDM – Mammography
Pattern Check (Used pattern)	-	-
Luminance Check	L'max / L'min > 170 L'max > 171 cd/m ² Lamb < Lmin / 1.5	L'max / L'min > 250 L'max > 250 cd/m ² Lamb < Lmin / 1.5
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 10 % of GSDF
Uniformity Check	-	-
Multi-monitor	-	-

Pattern Check

The Shades of RGB pattern displays 18 gradation levels for each of Red, Green, and Blue for checking. Monochrome monitors cannot run (display) this pattern even if it has been specified as a display pattern.

The Bi-Weekly setting is not available in RadiCS. Specify Weekly instead. The Visual Check settings are the same as those for Bi-Weekly.

Luminance Check

$Lamb < Lmin/1.5$ is added in accordance with AAPM On-line Report No. 03.

Grayscale Check

Added in accordance with AAPM On-line Report No. 03.

Uniformity Check

Added in accordance with AAPM On-line Report No. 03.

Sensors

All the measurement devices can be used in accordance with AAPM On-line Report No. 03.

Cautions

As the guideline does not contain any description of the acceptance test, the same settings as those for the consistency test (annually) are configured.

NYC-PDM**RadiCS Setup**

	Acceptance Test / Consistency Test [Annually]		
	NYC PDM – Hospitals	NYC PDM – Clinical sites	NYC PDM – Mammography
Pattern Check (Used pattern)	-	-	-
Luminance Check	L'max / L'min > 250 L'max > 350 cd/m ² Lamb < Lmin / 1.5	L'max / L'min > 250 L'max > 250 cd/m ² Lamb < Lmin / 1.5	L'max / L'min > 250 L'max > 420 cd/m ² Lamb < Lmin / 1.5
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 10 % of GSDF	Target error rate < 10 % of GSDF
Uniformity Check	Grayscale: 204, 26 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 30 %	Grayscale: 204, 26 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 30 %	Grayscale: 204, 26 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 30 %
Multi-monitor	-	-	-

	Consistency Test [Bi-Weekly]		
	NYC PDM – Hospitals	NYC PDM – Clinical sites	NYC PDM – Mammography
Pattern Check (Used pattern)	Black SMPTE Shades of RGB White	Black SMPTE Shades of RGB White	Black SMPTE Shades of RGB White
Luminance Check	-	-	-
Grayscale Check	-	-	-
Uniformity Check	-	-	-
Multi-monitor	-	-	-

	Consistency Test [Quarterly]		
	NYC PDM – Hospitals	NYC PDM – Clinical sites	NYC PDM – Mammography
Pattern Check (Used pattern)	-	-	-
Luminance Check	L'max / L'min > 250 L'max > 350 cd/m ² Lamb < Lmin / 1.5	L'max / L'min > 250 L'max > 250 cd/m ² Lamb < Lmin / 1.5	L'max / L'min > 250 L'max > 420 cd/m ² Lamb < Lmin / 1.5
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 10 % of GSDF	Target error rate < 10 % of GSDF
Uniformity Check	-	-	-
Multi-monitor	-	-	-

Pattern Check

The Shades of RGB pattern displays 18 gradation levels for each of Red, Green, and Blue for checking. Monochrome monitors cannot run (display) this pattern even if it has been specified as a display pattern.

The Bi-Weekly setting is not available in RadiCS. Specify Weekly instead. The Visual Check settings are the same as those for Bi-Weekly.

Luminance Check

$L_{amb} < L_{min}/1.5$ is added in accordance with AAPM On-line Report No. 03.

Grayscale Check

Added in accordance with AAPM On-line Report No. 03.

Uniformity Check

Each judgment condition includes an equality sign, but RadiCS does not.

Sensors

All the measurement devices can be used in accordance with AAPM On-line Report No. 03.

Cautions

The judgment of the luminance check has been added to each test. In addition, the judgment of the luminance ratio has been added to consistency tests (quarterly).

ONR 195240-20**RadiCS Setup**

	Acceptance Test			
	Category A	Category A Mammo	Category B	Category B Dentistry
Pattern Check (Used pattern)	TG18-OIQ TG18-UN80 TG18-UN10	TG18-OIQ TG18-UN80 TG18-UN10 TG18-MM1 TG18-MM2	TG18-OIQ TG18-UN80 TG18-UN10	TG18-OIQ TG18-UN80 TG18-UN10
Illuminance judgment	$\leq 50 \text{ lx}$	$\leq 50 \text{ lx}$	$\leq 100 \text{ lx}$	$\leq 100 \text{ lx}$
Luminance Check	$L'_{\max} / L'_{\min} > 100$ $L'_{\max} > 200 \text{ cd/m}^2$ $L_{\text{amb}} < L'_{\max} / 100$	$L'_{\max} / L'_{\min} > 250$ $L'_{\max} > 250 \text{ cd/m}^2$ $L_{\text{amb}} < L'_{\max} / 100$	$L'_{\max} / L'_{\min} > 40$ $L'_{\max} > 120 \text{ cd/m}^2$ $L_{\text{amb}} < L'_{\max} / 40$	$L'_{\max} / L'_{\min} > 40$ $L'_{\max} > 120 \text{ cd/m}^2$ $L_{\text{amb}} < L'_{\max} / 40$
Grayscale Check	-	-	-	-
Uniformity Check	Grayscale: 204, 26 ($L_{\text{corner}} - L_{\text{center}} / L_{\text{center}} \times 100 < 25 \%$)	Grayscale: 204, 26 ($L_{\text{corner}} - L_{\text{center}} / L_{\text{center}} \times 100 < 25 \%$)	Grayscale: 204, 26 ($L_{\text{corner}} - L_{\text{center}} / L_{\text{center}} \times 100 < 30 \%$)	Grayscale: 204, 26 ($L_{\text{corner}} - L_{\text{center}} / L_{\text{center}} \times 100 < 30 \%$)
Multi-monitor	$\Delta L'_{\max} < 20 \%$	$\Delta L'_{\max} < 10 \%$	$\Delta L'_{\max} < 20 \%$	$\Delta L'_{\max} < 20 \%$

	Consistency Test			
	Category A	Category A Mammo	Category B	Category B Dentistry
Pattern Check (Used pattern)	TG18-OIQ TG18-UN80 TG18-UN10	TG18-OIQ TG18-UN80 TG18-UN10 TG18-MM1 TG18-MM2	TG18-OIQ TG18-UN80 TG18-UN10	TG18-OIQ TG18-UN80 TG18-UN10
Illuminance judgment	$\leq 50 \text{ lx}$	$\leq 50 \text{ lx}$	$\leq 100 \text{ lx}$	-
Luminance Check	$L'_{\max} / L'_{\min} > 100$ $L'_{\max} > 200 \text{ cd/m}^2$ $L_{\text{amb}} < L'_{\max} / 100$ $\Delta L_{\text{amb}} < 30 \%$	$L'_{\max} / L'_{\min} > 250$ $L'_{\max} > 250 \text{ cd/m}^2$ $L_{\text{amb}} < L'_{\max} / 100$ $\Delta L_{\text{amb}} < 30 \%$	$L'_{\max} / L'_{\min} > 40$ $L'_{\max} > 120 \text{ cd/m}^2$ $L_{\text{amb}} < L'_{\max} / 40$ $\Delta L_{\text{amb}} < 30 \%$	-
Grayscale Check	-	-	-	-
Uniformity Check	Grayscale: 204, 26 ($L_{\text{corner}} - L_{\text{center}} / L_{\text{center}} \times 100 < 25 \%$)	Grayscale: 204, 26 ($L_{\text{corner}} - L_{\text{center}} / L_{\text{center}} \times 100 < 25 \%$)	Grayscale: 204, 26 ($L_{\text{corner}} - L_{\text{center}} / L_{\text{center}} \times 100 < 30 \%$)	-
Multi-monitor	$\Delta L'_{\max} < 20 \%$	$\Delta L'_{\max} < 10 \%$	$\Delta L'_{\max} < 20 \%$	-

ONR 195240-20: 2008 and RadiCS

Pattern Check

RadiCS prepares the patterns based on check results for respective compatible resolutions.

Luminance Check

L_{max} and L_{min} in ONR 195240-20, which include ambient luminance, are equivalent to L'_{max} and L'_{min} in RadiCS. L_{amb} indicates ambient luminance, the same value as "Ls" in ONR 195240-20. The equation is transformed by changing $L_{max}/L_s > 100$ (or 40) in ONR 195240-20 into $L_s < L_{max}/100$ (or 40). When a contact sensor is used in the monitor equipped with the illuminance sensor capable of measuring environmental illumination (see [5 Monitör Durumunu Kontrol Etme \[► 102\]](#)), the conversion from illuminance to brightness is automatically performed.

In RadiCS, as in accordance with the standard, no judgment will be made for Delta L_{amb} if the measurement value of the consistency test is 0.15 cd/m^2 or less and below the baseline value.

Uniformity Check

Luminance uniformity is determined from the ratio of difference in luminance between the center of the screen and a corner, with the center as the standard. ONR 195240-20 provides a method that uses the SMPTE pattern and another method that uses the TG18-UNL80 (or UNL10). RadiCS adopts the method that uses the TG18-UNL80 (or UNL10) pattern. It displays grayscale 204 and grayscale 26 windows (a square occupying 10 % of the total display area) in the center of the screen and corners, and measures the middle portion of the window.

All monitors compatible with RadiCS are LCD, therefore, LCD values (25 % and 30 %) are used as the judgment value. For this reason, CRT monitors are not supported.

RadiCS specifies $(L_{corner} - L_{center})/L_{center} \times 100 < 25 \%$ (or 30 %), but this denotes $\pm 25 \%$ (or $\pm 30 \%$), and does not include an equals sign.

Sensors

For acceptance tests, ONR 195240-20 defines the use of measurement devices conforming to class B or higher (DIN 5032-7) and those do not block ambient light. To perform acceptance tests using RadiCS, only non-contact type measurement devices can be used. EIZO sensors can also be used for consistency tests.

Multi-monitor

ONR 195240-20 has multi-monitor judgment. If necessary, make settings as indicated in the table above. ONR 195240-20 includes an equality sign but RadiCS does not.

Cautions

Category A Mammo requires a minimum resolution of 2000×2500 for monitors used for mammography, however, RadiCS does not perform this judgment.

DIN 6868-157**RadiCS Setup**

	Acceptance Test				
	I. Mammography	II. Mammographic stereotaxy	III. Projection radiography	IV. Fluoroscopy, all applications	V. Computed tomography
Pattern Check (Used pattern)	TG18-OIQ TG18-UN80 TG18-UN10 TG18-MP TG18-LPH (89,50,10) TG18-LPV (89,50,10)	TG18-OIQ TG18-UN80 TG18-UN10 TG18-MP			
Luminance Check	L'max > 250 cd/m ² L'max / L'min > 250	L'max > 200 cd/m ² L'max / L'min > 100	L'max > 250 cd/m ² L'max / L'min > 250	L'max > 150 cd/m ² L'max / L'min > 100	
	Lamb < Lmin / 0.1 ^{*1}				
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 15 % of GSDF	Target error rate < 10 % of GSDF	Target error rate < 15 % of GSDF	
Uniformity Check	Grayscale: 204, 26 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 25 %				
Multi-monitor ^{*2}	Grayscale: 26 < 10 %	Grayscale: 26 < 20 %			
Resolution	≥2048 x ≥2048	≥1024 x ≥1024	≥1600 x ≥1200	≥1024 x ≥1024	≥1024 x ≥1024

	Acceptance Test		
	VI. Digital volume tomography (dental) etc. in RK 5	VII. Intraoral X-ray diagnostics (dental) etc. in RK 6	VIII. Viewing
	VI. Dental X-ray equipment etc. in RK 5 (five-year interval) ^{*3}		
Pattern Check (Used pattern)	TG18-OIQ TG18-UN80 TG18-UN10		-
Luminance Check	L'max > 200 cd/m ² L'max / L'min > 100	L'max > 300 cd/m ² L'max / L'min > 100	-
	$\text{Lamb} < \text{Lmin} / 0.1^{*1}$		-
Grayscale Check	-	-	-
Uniformity Check	Grayscale: 204, 26 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 30 %		-
Multi-monitor ^{*2}	Grayscale: 26 < 30 %		-

	Acceptance Test		
	VI. Digital volume tomography (dental) etc. in RK 5	VII. Intraoral X-ray diagnostics (dental) etc. in RK 6	VIII. Viewing
	VI. Dental X-ray equipment etc. in RK 5 (five-year interval)* ³		
Resolution	≥1024 x ≥768		-

*1 $L'_{min} > 1.1L_{amb} \leftrightarrow L_{min} + L_{amb} > 1.1L_{amb} \leftrightarrow L_{min} > 0.1L_{amb} \leftrightarrow L_{amb} < L_{min}/0.1$

*2 $(L_{high} - L_{low}) / (L_{high} + L_{low}) \times 200$

*3 The interval of the annual measuring tests can be extended to five years on the condition that the specified requirements are satisfied.

	Consistency Test				
	I. Mammography	II. Mammographic stereotaxy	III. Projection radiography	IV. Fluoroscopy, all applications	V. Computed tomography
Pattern Check (Used pattern)	TG18-OIQ TG18-UN80				
Luminance Check	L'max > 250 cd/m ²	L'max > 200 cd/m ²	L'max > 250 cd/m ²	L'max > 150 cd/m ²	
	L'max / L'min > 250	L'max / L'min > 100	L'max / L'min > 250	L'max / L'min > 100	
	Lamb < Lmin / 0.1 ^{*1}				
	ΔL'max < 30 %				
	ΔL'min < 30 %				
	ΔLamb ≤ 30 % ^{*3}				
Grayscale Check ^{*3}	Target error rate < 10 % of GSDF	Target error rate < 15 % of GSDF	Target error rate < 10 % of GSDF	Target error rate < 15 % of GSDF	
Uniformity Check	-	-	-	-	-
Multi-monitor ^{*2, 3}	Grayscale: 26 < 10 %	Grayscale: 26 < 20 %			
Resolution	≥2048 x ≥2048	≥1024 x ≥1024	≥1600 x ≥1200 ^{*4}	≥1024 x ≥1024	≥1024 x ≥1024

	Consistency Test		
	VI. Digital volume tomography (dental) etc. in RK 5	VII. Intraoral X-ray diagnostics (dental) etc. in RK 6	VIII. Viewing
	VI. Dental X-ray equipment etc. in RK 5 (five-year interval) ⁵		
Pattern Check (Used pattern)	TG18-OIQ TG18-UN80		TG18-OIQ
Luminance Check	L'max > 200 cd/m ²	L'max > 300 cd/m ²	-
	L'max / L'min > 100	L'max / L'min > 100	
	Lamb < Lmin / 0.1 ^{*1} $\Delta L'max < 30\%$ $\Delta L'min < 30\%$		-
	-	-	-
Grayscale Check	-	-	-
Uniformity Check	-	-	-
Multi-monitor ^{*2}	-	-	-
Resolution	$\geq 1024 \times \geq 768$		-

*1 $L'min \geq 1.1Lamb \leftrightarrow Lmin + Lamb \geq 1.1Lamb \leftrightarrow Lmin \geq 0.1Lamb \leftrightarrow Lamb \leq Lmin/0.1$

*2 $(L_{high} - L_{low}) / (L_{high} + L_{low}) \times 200$

*3 If Room Category "RK3" is selected, it will be excepted from judgment. If the luminance satisfies $\Delta L'min < 30\%$, $\Delta Lamb < 30\%$ does not display or provide judgment.

*4 $\geq 1024 \times \geq 1024$ can be used until December 31, 2024 as per transition measures.

*5 The interval of the annual measuring tests can be extended to five years on the condition that the specified requirements are satisfied.

Correlation between DIN 6868-157 and RadiCS

Test requirements

To create a test result report in RadiCS, it is necessary to check and enter information of the requirements before executing the test.

- Check that the image display system has adequate ability and has been installed and configured correctly.
(E.g., the system is for medical use, the grayscale characteristics of the image display device are GSDF, and the system has been configured and installed correctly according to the specifications.)
- Check that the specifications of the measurement device and software to be used in the test are appropriate.
(E.g., using the measurement device of DIN 5032-7 class B or higher, acceptance test, selecting and securing the reference clinical image^{*1}, appropriate resolution of the test image^{*2}, ensuring reliability of the testing software^{*2})
- Check that the environment where the test is executed has been set up.
(E.g., turn on the power of the monitor in advance, clean the display, stabilize the ambient light, and prevent reflection.)

As DIN 6868-157 specifies not only selection of the body parts and capture methods but also illuminance that should be selected depending on the actual work and locations, so it is necessary to select the environmental illumination^{*3}. RK that can be selected differs depending on the selected body part and capture method.

Room category	Location (Work)	Illuminance (lx)
RK1	Diagnostics room	≤50
RK2	Examination rooms with immediate diagnostics	≤100
RK3	Rooms to carry out examinations	≤500
RK4	Viewing and treatment rooms	≤1000
RK5	Dental diagnostic workstation	≤100
RK6	Dental treatment room	≤1000

- *1 An appropriate clinical image should be selected as reference clinical image and viewed with optimum parameters. Before running RadiCS, check the quality of the image secured by the responsible operator on the application software (viewer, etc.) to be actually used for displaying the image. On the reference clinical image confirmation dialog, enter the image identification, parameters to be displayed, name of the responsible operator, and other necessary information. Enter the judgment result when performing pattern check.
- *2 RadiCS displays the test image in the same resolution as that of the monitor, so each pixel of the test image corresponds to that of the monitor. As displayed image is not corrected by the software, it is possible to evaluate the monitor characteristics correctly even in measurement of grayscale characteristics such as GSDF.
- *3 It may be necessary for the environmental illumination to be set appropriately in order to pass the test.

Pattern Check

RadiCS determines the properties to be verified and independently prepares patterns for each compatible resolution applied.

As for checking the reference clinical image, the items to be checked are displayed but the image is not displayed. As the check here is only for recording the history of check results, you need to judge with the results you checked before execution.

Although the TG18-MP pattern has been created as a pattern of 10 bits or more enabling identification of both 8 bits and 10 bits resolutions, RadiCS creates and displays it as an 8-bit pattern. An 8-bit pattern is enough to check the judgment criteria of the test items.

Luminance Check

In case of DIN 6868-157, luminance of ambient light should be included in the test. When a contact sensor is used in the monitor equipped with the illuminance sensor capable of measuring environmental illumination (see [5 Monitör Durumunu Kontrol Etme](#) ► 102), the conversion from illuminance to brightness is automatically performed.

Deviation from the reference value includes an equality sign in the standard but not in RadiCS.

Therefore $L'_{min} \geq 1.1L_{amb}$ does not include an equality sign in RadiCS.

In RadiCS, as in accordance with the standard, no judgment will be made for Delta Lamb if the measurement value of the consistency test is 0.15 cd/m² or less and below the baseline value.

Grayscale Check

GSDF checking includes an equality sign in the standard but not in RadiCS.

Uniformity Check

In DIN 6868-157, luminance uniformity is measured at five points for less than 23 inch and nine points for 23 inch or larger, which will be selected automatically.

If a contact type sensor is used, luminance of ambient light is not included.

" $(L_{max} - L_{min}) / (L_{max} + L_{min}) \times 200$ " shown in RadiCS is the same as " $200 \times (L_{highest} - L_{lowest}) / (L_{highest} + L_{lowest})$ " in the standard.

Sensors

DIN 6868-157 requires a luminance meter class B or higher (DIN 5032-7) for acceptance tests and measurement devices that does not block environmental light. If measuring grayscale by bringing a measurement device in contact with the monitor, use a measurement device that, in accordance with the measurement devices' User's Manual, can be brought in contact with the monitor.

EIZO sensors are available for consistency tests. DIN 6868-157 requires creation of a reference value for consistency test to include reflected luminance caused by ambient light and accepts the use of contact type sensor.

If any measurement device or measurement method different from that is used in the acceptance test is used, it is recommended to make a correlation with the measurement device used in the acceptance test before deciding the reference value.

Multi-monitor

DIN 6868-157 includes a determination for multiple monitors, but by default RadiCS is set not to make such a determination. Enter the settings as necessary (see [RadiCS Setup](#) [▶ 212]).

"(Lhigh - Llow)/(Lhigh + Llow) x 200" shown in RadiCS is the same as "200 x (Lhighest - Llowest)/(Lhighest + Llowest)" in the standard.

Resolution

The available monitor resolution is determined in the standard depending on body part / capture method. RadiCS has set restrictions in the control criteria to be selected for body parts / capture methods in accordance with the standard.

	I. Mammography	II. Mammographic stereotaxy	III. Projection radiography	IV. Fluoroscopy, all applications / V. Computed tomography	VI. Digital volume tomography (dental) etc. in RK 5/ VI. Dental X-ray equipment etc. in RK 5 (five-year interval) / VII. Intraoral X-ray diagnostics (dental) etc. in RK 6
Resolution	≥2048 x ≥2048	≥1024 x ≥1024	≥1600 x ≥1200	≥1024 x ≥1024	≥1024 x ≥768

DIN V 6868-57**RadiCS Setup**

	Acceptance Test	
	Category A	Category B
Pattern Check (Used pattern)	Test pattern 1 Test pattern 2 Test pattern 3	Test pattern 1 Test pattern 2 Test pattern 3
Luminance Check	$L'max / L'min > 100$ $Lamb < L'max / 100$	$L'max / L'min > 40$ $Lamb < L'max / 40$
Grayscale Check	-	-
Uniformity Check	Grayscale: 128 ^{*1}	Grayscale: 128 ^{*2}

	Consistency Test	
	Category A	Category B
Pattern Check (Used pattern)	Test pattern 1 Test pattern 2 Test pattern 3	Test pattern 1 Test pattern 2 Test pattern 3
Luminance Check	$L'max / L'min > 100$ $Lamb < L'max / 100$	$L'max / L'min > 40$ $Lamb < L'max / 40$
Grayscale Check	-	-
Uniformity Check	Grayscale: 128 ^{*1}	Grayscale: 128 ^{*2}

*1 $(L_{corner} - L_{center}) / L_{center} \times 100 < 15 \%$

*2 $(L_{corner} - L_{center}) / L_{center} \times 100 < 20 \%$

Correlation between DIN V 6868-57 and RadiCS**Pattern Check**

A test pattern given in DIN V 6868-57 cannot be applied to a monitor whose screen aspect ratio is not 1:1 without modification, since DIN V 6868-57 (or the test pattern) uses an aspect ratio of 1:1. Therefore, RadiCS checks a monitor being tested, and determines and generates an appropriate test pattern for each resolution supported by the monitor.

- Test pattern 1
Equivalent to Bild 3 pattern. The pattern is scaled in accordance with the screen resolution.
- Test pattern 2
Equivalent to Bild 2 pattern. The pattern is scaled in accordance with the screen resolution.
- Test pattern 3
Equivalent to Bild 5 pattern. The pattern is scaled in accordance with the screen resolution.

Luminance Check

L_{max} and L_{min} used in DIN V 6868-57 include ambient luminance and are the same as L'_{max} and L'_{min} in RadiCS. L_{amb} stands for the ambient luminance and refers to the same value as "L_s" of DIN V 6868-57. L_{max}/L_s > 100 (or 40) have been L_s < L_{max}/100 (or 40).

L'_{max}/L'_{min} stands for a contrast ratio. DIN V 6868-57 includes an equality sign like L_{max}/L_{min} ≥ 100 (or 40) but RadiCS does not.

DIN V 6868-57 defines L'_{max} and L'_{min} by measuring the test pattern 2 square with white (grayscale: 255) and black (grayscale: 0). RadiCS displays 10 % of a display area in the middle and measures luminance by changing the grayscale 0 to 255. By doing so, the exact contrast ratio can be acquired.

Uniformity Check

The Uniformity Check judges the uniformity of the ratio between the screen corner and the center of the screen as a standard. DIN V 6868-57 has no particular standard regarding measuring points. It also displays 10 % display area of the window at grayscale 128 in the middle of the screen and in the corner of the screen and measures the center of the window.

The basic judgment value (15 % or 20 %) is the same as LCD monitors since RadiForce series monitors are recommended for RadiCS.

RadiCS specifies (L_{corner}-L_{center})/L_{center} x 100 < 15 % (or 20 %), but this denotes ±15 % (or ±20 %), and does not include an equals sign.

Sensors

DIN V 6868-57 requires a luminance meter class B or higher (DIN 5032-7) for acceptance tests and measurement devices that does not block environmental light.

DIN V 6868-57 allows noncontact sensors only to measure Category B reference value for consistency tests. EIZO sensors are available for consistency tests.

DMG QC Manual

RadiCS Setup

	Acceptance Test	Consistency Test
Pattern Check (Used pattern)	TG18-QC TG18-UN80	Black TG18-QC TG18-UN80
Luminance Check	$L'_{\max} / L'_{\min} > 250$	$L'_{\max} / L'_{\min} > 250$ $\Delta L'_{\max} < 10 \%$
Grayscale Check	Target error rate $< 15 \%$ of GSDF	Target error rate $< 15 \%$ of GSDF
Uniformity Check	Grayscale: 204 $(L_{\max} - L_{\min}) / (L_{\max} + L_{\min}) \times 200 < 30 \%$	-
Multi-monitor	$\Delta L'_{\max} < 10 \%$ between multiple monitors	$\Delta L'_{\max} < 10 \%$ between multiple monitors

Correlation between DMG QC Manual and RadiCS

Pattern Check

RadiCS determines necessary test patterns based on the inspection results and generates its own patterns corresponding to the resolution of the monitor.

- TG18-QC
Equivalent to the pattern with the same name in the standard. However, RadiCS-specific scaling is performed in accordance with the monitor resolution.
- TG18-UN80
A pattern solidly filled with white of grayscale 204. The pattern with the same name in the JESRA has a square frame, but RadiCS does not have it because it is not necessary for the visual inspection.

Luminance Check

In DMG QCM, the luminance measurement does not include the ambient luminance. In RadiCS, an apostrophe (') in the $L'_{\max}$, for example, indicates that it includes the ambient luminance. However, entering the ambient luminance value as 0 cd/m^2 can effectively exclude the ambient luminance from the luminance measurement.

Note that none of inequalities used in RadiCS includes an equality sign although every judgment condition in DMG QCM includes it.

The $L_{\max}$ value in the calibration setup is provided as the default for the baseline value of $\Delta L'_{\max}$.

Grayscale Check

In DMG QCM, the luminance measurement does not include the ambient luminance. In RadiCS, an apostrophe (') in the $L'_{\max}$, for example, indicates that it includes the ambient luminance. However, entering the ambient luminance value as 0 cd/m^2 can effectively exclude the ambient luminance from the luminance measurement.

The calculation method for this item is the same as the one for $\kappa\delta$. RadiCS describes the specification of the grayscale as Target Error Rate $< 15 \%$ (or 30%) of GSDF. Note that none of inequalities used in RadiCS includes an equality sign.

This specification is provided as the judgment condition for DICOM Part 14 GSDF, so there is no meaning to use this specification for other display functions. The number of measuring points is fixed to 18 points and this value cannot be changed. (The number of data points will be 17 because the result is presented as $(JND_{n+1} - JND_n)/2$.)

Uniformity Check

Although DMG QCM includes an equality sign, each judgment condition in RadiCS does not.

The DMG QCM specifies that the luminance is measured using the TG18-UN80 patterns. On the other hand, RadiCS displays two windows (grayscale: 204) with the size of 10 % of the whole display area at the center and a corner of the screen. It then measures the luminance at the center of both windows.

Sensors

DMG QCM permits the use of both noncontact and contact type measurement devices. In RadiCS, the noncontact measurement device measures the monitor without shutting off the environment light, so use the device in a dark room or use a cylinder to shut off environment light. Any sensors can be used to perform both the acceptance tests and the consistency tests.

Multi-monitor

DMG QCM has multi-monitor judgment. DMG QCM includes an equality sign but RadiCS does not.

EUREF**RadiCS Setup**

	Acceptance Test	
	Primary	Secondary
Pattern Check (Used pattern)	TG18-QC TG18-LPH (89, 50, 10) TG18-LPV (89, 50, 10)	TG18-QC TG18-LPH (89, 50, 10) TG18-LPV (89, 50, 10)
Luminance Check	$L'_{\max} / L'_{\min} > 250$	$L'_{\max} / L'_{\min} > 100$
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 20 % of GSDF
Uniformity Check	Grayscale: 26 $(L_{\max} - L_{\min}) / L_{\text{center}} \times 100 < 30 \%$	Grayscale: 26 $(L_{\max} - L_{\min}) / L_{\text{center}} \times 100 < 30 \%$
	Grayscale: 204 $(L_{\max} - L_{\min}) / L_{\text{center}} \times 100 < 15 \%$	Grayscale: 204 $(L_{\max} - L_{\min}) / L_{\text{center}} \times 100 < 15 \%$
Multi-monitor	$\Delta L'_{\max} < 5 \%$ between multiple monitors	$\Delta L'_{\max} < 5 \%$ between multiple monitors

	Consistency Test	
	Primary	Secondary
Pattern Check (Used pattern)	TG18-QC TG18-LPH (89, 50, 10) TG18-LPV (89, 50, 10)	TG18-QC TG18-LPH (89, 50, 10) TG18-LPV (89, 50, 10)
Luminance Check	$L'_{\max} / L'_{\min} > 250$	$L'_{\max} / L'_{\min} > 100$
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 20 % of GSDF
Uniformity Check	Grayscale: 26 $(L_{\max} - L_{\min}) / L_{\text{center}} \times 100 < 30 \%$	Grayscale: 26 $(L_{\max} - L_{\min}) / L_{\text{center}} \times 100 < 30 \%$
	Grayscale: 204 $(L_{\max} - L_{\min}) / L_{\text{center}} \times 100 < 15 \%$	Grayscale: 204 $(L_{\max} - L_{\min}) / L_{\text{center}} \times 100 < 15 \%$
Multi-monitor	$\Delta L'_{\max} < 5 \%$ between multiple monitors	$\Delta L'_{\max} < 5 \%$ between multiple monitors

Correlation between EUREF and RadiCS**Pattern Check**

The patterns used for EUREF are the same as those used for AAPM. RadiCS determines the properties to be verified and independently prepares appropriate patterns for each resolution.

- TG18-QC
This is scaled to match the resolution.
- TG18-LPH (89, 50, 10)
This is scaled to match the resolution.
- TG18-LPV (89, 50, 10)
This is scaled to match the resolution.

Luminance Check

Maximum luminance and luminance ratio specified in the standard correspond to L'max and L'max/L'min used in RadiCS. The patterns TG18-LN12-01 and TG18-LN12-18 are recommended for luminance measurements, but RadiCS measures the luminance by displaying a window equivalent to 10 % of the display area in the center of the screen and changing its grayscale level to 0 and 255. This provides a more accurate measurement. EUREF includes an equality sign but RadiCS does not.

Grayscale Check

The GSDF determination conditions correspond to those specified in EUREF. EUREF recommends using patterns TG18-LN12-01 to TG18-LN12-18 for measurements, but RadiCS measures the luminance by displaying a window equivalent to 10 % of the display area in the center of the screen and changing the grayscale level corresponding to the specified pattern from 0 to 255. This provides a more accurate measurement. EUREF includes an equality sign but RadiCS does not.

Uniformity Check

EUREF recommends using the TG18-UNL10 and TG18-UNL80 patterns, but since they have an aspect ratio of 1:1 they cannot be used directly. Instead, RadiCS displays grayscale 204 and grayscale 26 windows equivalent to 10 % of the display area in the center of the screen and in the corners, and measures the center portion of each window.

In Supplements: 2013, the judgment standard for LCDs to satisfy in relation to grayscale 204 has been tightened from 30 % to 15 % (30 % for CRTs). RadiCS monitors satisfy the standard applicable to LCDs.

Sensors

EUREF recommends the use of a telescopic luminance meter. EIZO sensors may also be used to perform measurements.

Multi-monitor

EUREF includes a determination for multiple monitors, but by default RadiCS is set not to make such a determination. If necessary, make settings as indicated in the table above. EUREF includes an equality sign but RadiCS does not.

Cautions

For primary use, an illuminance meter must be used to ensure that the ambient light level is less than 10 lux. RadiCS does not make illuminance-based judgment.

RadiForce series monitors are considered to sufficiently satisfy requirements regarding geometrical distortion, so this item is omitted.

IPEM

RadiCS Setup

	Acceptance Test	Consistency Test
Pattern Check (Used pattern)	TG18-QC	TG18-QC
Luminance Check	$L'max / L'min > 250$ $\Delta L'max < 20 \%$	$L'max / L'min > 250$ $\Delta L'max < 20 \%$
Grayscale Check	Target error rate $< 10 \%$ of GSDF	Target error rate $< 10 \%$ of GSDF
Uniformity Check	Grayscale: 128 $(Lmax-Lmin) / (Lmax+Lmin) \times 200 < 30 \%$	Grayscale: 128 $(Lmax-Lmin) / (Lmax+Lmin) \times 200 < 30 \%$
Multi-monitor	$\Delta L'max < 30 \%$ between multiple monitors $\Delta L'min < 30 \%$ between multiple monitors	$\Delta L'max < 30 \%$ between multiple monitors $\Delta L'min < 30 \%$ between multiple monitors

Correlation between IPEM and RadiCS

Pattern Check

The patterns used for IPEM are the same as those used for AAPM. RadiCS determines the properties to be verified and independently prepares appropriate patterns for each resolution.

- TG18-QC

This is scaled to match the resolution.

Luminance Check

Maximum luminance and luminance ratio specified in IPEM correspond to $L'max$ and $L'max/L'min$ used in RadiCS. The patterns TG18-QC and SMPTE are recommended for luminance measurements, but RadiCS measures the luminance by displaying a window equivalent to 10 % of the display area in the center of the screen and changing its grayscale level to 0 and 255. This provides a more accurate measurement. IPEM makes $\Delta Lmin \leq 25 \%$ judgment, but RadiCS does not. Make the settings as necessary although the standard name will be "Custom". IPEM includes an equality sign but RadiCS does not.

Uniformity Check

IPEM recommends using TG18-QC or SMPTE patterns, but these patterns are not suitable for measuring 50 % grayscale uniformity. Instead, RadiCS displays grayscale 128 windows equivalent to 10 % of the display area in the center of the screen and in the corners, and measures the center portion of each window. IPEM includes an equality sign but RadiCS does not.

Sensors

Use of a measurement device that complies with the CIE standard photopic spectral response and has a calibration traceable to an appropriate primary standard is recommended. RadiCS supports use of all compliant sensors.

Multi-monitor

IPEM includes a determination for multiple monitors, but by default RadiCS is set not to make such a determination. If necessary, make settings as indicated in the table above. IPEM includes an equality sign but RadiCS does not.

Cautions

An illuminance meter must be used to ensure that the ambient light level is less than 15 lux. RadiCS does not make illuminance-based judgment.

JESRA TR-0049**RadiCS Setup**

	Acceptance Test				
	Category I-A	Category I-B	Category II (for diagnosis plus)	Category II (for diagnosis)	Category II (for reference)
Pattern Check (Used pattern)	TG18-OIQ TG18-MP TG18-UN80 TG18-UN10 JESRA Clinical Image	TG18-OIQ TG18-MP TG18-UN80 TG18-UN10 JESRA Clinical Image	TG18-OIQ TG18-MP TG18-UN80 JESRA Clinical Image	TG18-OIQ TG18-MP TG18-UN80 JESRA Clinical Image	TG18-OIQ TG18-MP TG18-UN80 JESRA Clinical Image
Luminance Check	L'max > 450 cd/m ² L'max / L'min > 350 $\Delta L'max < \pm 10 \%$ $Lamb < Lmin / 0.67$	L'max > 350 cd/m ² L'max / L'min > 250 $\Delta L'max < \pm 10 \%$ $Lamb < Lmin / 0.67$	L'max > 170 cd/m ² L'max / L'min > 250 $\Delta L'max < \pm 10 \%$	L'max > 150 cd/m ² L'max / L'min > 100 $\Delta L'max < \pm 10 \%$	L'max > 150 cd/m ² L'max / L'min > 100 $\Delta L'max < \pm 10 \%$
Grayscale Check	Target error rate < 10 % of GSDF Grayscale chromaticity $\Delta u'v' < 0.010$ (5.00 cd/m ² or more)	Target error rate < 10 % of GSDF Grayscale chromaticity $\Delta u'v' < 0.010$ (5.00 cd/m ² or more)	Target error rate < 15 % of GSDF Grayscale chromaticity $\Delta u'v' < 0.015$ (5.00 cd/m ² or more)	Target error rate < 20 % of GSDF Grayscale chromaticity $\Delta u'v' < 0.015$ (5.00 cd/m ² or more)	Target error rate < 20 % of GSDF
Uniformity Check	Grayscale 204 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 20 % Grayscale 204 $\Delta u'v' < 0.010$	Grayscale 204 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 20 % Grayscale 204 $\Delta u'v' < 0.010$	Grayscale 204 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 30 % Grayscale 204 $\Delta u'v' < 0.010$	Grayscale 204 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 30 % Grayscale 204 $\Delta u'v' < 0.015$	Grayscale 204 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 30 %
Multi-monitor	$\Delta L'max < 10 \%$ $\Delta u'v' < 0.010$	$\Delta L'max < 10 \%$ $\Delta u'v' < 0.010$	$\Delta L'max < 10 \%$ $\Delta u'v' < 0.010$	$\Delta L'max < 20 \%$ $\Delta u'v' < 0.015$	$\Delta L'max < 20 \%$
Resolution	$\geq 2048 \times \geq 2048$	$\geq 1024 \times \geq 1024$	$\geq 1000 \times \geq 1000$	-	-

	Consistency Test				
	Category I-A	Category I-B	Category II (for diagnosis plus)	Category II (for diagnosis)	Category II (for reference)
Pattern Check (Used pattern)	TG18-OIQ TG18-MP TG18-UN80 TG18-UN10 JESRA Clinical Image	TG18-OIQ TG18-MP TG18-UN80 TG18-UN10 JESRA Clinical Image	TG18-OIQ TG18-MP TG18-UN80 JESRA Clinical Image	TG18-OIQ TG18-MP TG18-UN80 JESRA Clinical Image	TG18-OIQ TG18-MP TG18-UN80 JESRA Clinical Image
Luminance Check	$L'_{max} > 450 \text{ cd/m}^2$ $L'_{max} / L'_{min} > 350$ $\Delta L'_{max} < \pm 10 \%$ $L_{amb} / L_{min} < 0.67$	$L'_{max} > 350 \text{ cd/m}^2$ $L'_{max} / L'_{min} > 250$ $\Delta L'_{max} < \pm 10 \%$ $L_{amb} / L_{min} < 0.67$	$L'_{max} > 170 \text{ cd/m}^2$ $L'_{max} / L'_{min} > 250$ $\Delta L'_{max} < \pm 10 \%$	$L'_{max} > 150 \text{ cd/m}^2$ $L'_{max} / L'_{min} > 100$ $\Delta L'_{max} < \pm 10 \%$	$L'_{max} > 150 \text{ cd/m}^2$ $L'_{max} / L'_{min} > 100$ $\Delta L'_{max} < \pm 10 \%$
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 10 % of GSDF	Target error rate < 15 % of GSDF	Target error rate < 20 % of GSDF	Target error rate < 20 % of GSDF
Uniformity Check	-	-	-	-	-
Multi-monitor	$\Delta L'_{max} < 10 \%$	$\Delta L'_{max} < 10 \%$	$\Delta L'_{max} < 10 \%$	$\Delta L'_{max} < 20 \%$	$\Delta L'_{max} < 20 \%$
Resolution	$\geq 2048 \times \geq 2048$	$\geq 1024 \times \geq 1024$	$\geq 1000 \times \geq 1000$	-	-

Correlation between JESRA TR-0049 and RadiCS

Pattern Check

RadiCS prepares the patterns based on check results for respective compatible resolutions. In RadiCS, the test pattern is labeled as "TG18-OIQ," but it is identical in specifications to the "OIQ" test pattern.

Luminance Check

The standard includes an equality sign in each judgment condition but RadiCS does not include an equality sign.

The " L_{amb} / L'_{min} (a) relationship < 0.6" equation has been changed to " $L_{amb} < L_{min} / 0.67$ " to determine the ambient luminance.

Grayscale Check

The standard includes an equality sign in each judgment condition but RadiCS does not include an equality sign.

In RadiCS, "target error rate < 10 (15 %, 20 %) of GSDF" indicates a contrast response test, which measures 18 points. Measured values of less than 5.00 cd/m² are not used to determine "Grayscale chromaticity $\Delta u'v'$ ".

Uniformity Check

The standard includes an equality sign in each judgment condition but RadiCS does not include an equality sign.

In JESRA TR-0049, measurements are performed while displaying the TG18-UNL80 pattern on the full screen. In RadiCS, window patterns (same as the TG18-UN80 specifications), each of which is 10 % of the display area in 204 gradations, are sequentially displayed in the center or corner of the screen, which enables an easy-to-perform measurement. In RadiCS, it is expressed $(L_{max} - L_{min}) / (L_{max} + L_{min}) \times 200$.

Sensors

Noncontact (telescopic) and contact measurement devices are available in JESRA TR-0049.

Multi-monitor

The standard includes multi-monitor judgment and includes an equality sign, but RadiCS does not include an equality sign.

JESRA X-0093**RadiCS Setup**

An apostrophe (') in L'max and L'min indicates that it includes the ambient luminance. However, using a measurement method that does not include the ambient luminance or by entering the ambient luminance value as "0 cd/m²", judgment can exclude the ambient luminance from the luminance measurement.

In RadiCS, each condition does not include this symbol; however, this fact has no real influence because judgment is performed using a lower value than the fourth decimal place.

	Acceptance Test		
	Grade 1A	Grade 1B	Grade 2
Pattern Check (Used pattern)	TG18-QC TG18-UN80 JESRA Clinical Image	TG18-QC TG18-UN80 JESRA Clinical Image	TG18-QC TG18-UN80 JESRA Clinical Image
Luminance Check	L'max / L'min > 250 L'max > 350 cd/m ²	L'max / L'min > 250 L'max > 170 cd/m ²	L'max / L'min > 100 L'max > 100 cd/m ²
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 15 % of GSDF	Target error rate < 30 % of GSDF
Uniformity Check	Grayscale: 204 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 30 % Grayscale: 204 $\Delta u'v' < 0.010$	Grayscale: 204 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 30 % Grayscale: 204 $\Delta u'v' < 0.010$	Grayscale: 204 (Lmax-Lmin) / (Lmax+Lmin) x 200 < 30 %
Multi-monitor	$\Delta L'max < 10 \%$ between multiple monitors Grayscale: 204 Mean value between multiple monitors $\Delta u'v' < 0.010$	$\Delta L'max < 10 \%$ between multiple monitors Grayscale: 204 Mean value between multiple monitors $\Delta u'v' < 0.010$	$\Delta L'max < 10 \%$ between multiple monitors

	Consistency Test		
	Grade 1A	Grade 1B	Grade 2
Pattern Check (Used pattern)	TG18-QC TG18-UN80 JESRA Clinical Image	TG18-QC TG18-UN80 JESRA Clinical Image	TG18-QC TG18-UN80 JESRA Clinical Image
Luminance Check	L'max / L'min > 250 L'max > 350 cd/m ² $\Delta L'max < 10 \%$	L'max / L'min > 250 L'max > 170 cd/m ² $\Delta L'max < 10 \%$	L'max / L'min > 100 L'max > 100 cd/m ² $\Delta L'max < 10 \%$
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 15 % of GSDF	Target error rate < 30 % of GSDF
Uniformity Check	-	-	-
Multi-monitor	$\Delta L'max < 10 \%$ between multiple monitors	$\Delta L'max < 10 \%$ between multiple monitors	$\Delta L'max < 10 \%$ between multiple monitors

Correlation between JESRA X-0093 and RadiCS**Pattern Check**

The guideline introduces test patterns for conducting a test, but it does not cover all medical monitors' resolutions. RadiCS provides the appropriate test patterns, taking into account the check contents shown in the guideline.

Luminance Check

The ambient change ratio between the baseline value and the measured value is indicated by " ΔL_{max} ". The default baseline value is set to the L_{max} value in the Calibration Settings.

Grayscale Check

The maximum error rate of contrast response, " $\kappa\delta$ ", is indicated by "target error rate < 10 % (15 %, 30 %) of GSDF".

Uniformity Check

In JESRA X-0093, measurements are performed while displaying the TG18-UN80 pattern on the full screen. In RadiCS, window patterns (same as the TG18-UN80 specifications), each of which is 10 % of the display area in 204 gradations, are sequentially displayed in the center or corner of the screen, which enables an easy-to-perform measurement. In RadiCS, the luminance uniformity is indicated by " $(L_{\text{max}} - L_{\text{min}})/(L_{\text{max}} + L_{\text{min}}) \times 200$ ".

Sensors

JESRA X-0093 provides use of both the non-contact type (telescopic) and contact type measurement devices; therefore, all the compatible sensors can be used.

The non-contact type measurement device performs measurements including the ambient luminance. When you do not want to include the ambient luminance, perform measurements in a dark room or shut down the environmental light using a circular cylinder, etc.

Multi-monitor

The differential ratio of the maximum luminance between medical monitors is indicated by " ΔL_{max} ".

QS-RL**RadiCS Setup**

	Acceptance Test		
	Category A	Category B	Category A Mammo
Pattern Check (Used pattern)	Test pattern 1 Test pattern 2 Test pattern 3	Test pattern 1 Test pattern 2 Test pattern 3	Test pattern 1 Test pattern 2 Test pattern 3
Luminance Check	$L'max / L'min > 100$ $L'max > 200 \text{ cd/m}^2$ $Lamb < L'max / 100$	$L'max / L'min > 40$ $L'max > 120 \text{ cd/m}^2$ $Lamb < L'max / 40$	$L'max / L'min > 250$ $L'max > 250 \text{ cd/m}^2$ $L'min > 1.0 \text{ cd/m}^2$ $Lamb < L'max / 100$
Grayscale Check	-	-	-
Uniformity Check	Grayscale: 128 (Lcorner-Lcenter) / Lcenter x 100 < 15 %	Grayscale: 128 (Lcorner-Lcenter) / Lcenter x 100 < 20 %	Grayscale: 128 (Lcorner-Lcenter) / Lcenter x 100 < 15 %
Multi-monitor	-	-	$\Delta L'max < 10 \%$ between multiple monitors $\Delta(L'max / L'min) < 10 \%$ between multiple monitors

	Consistency Test		
	Category A	Category B	Category A Mammo
Pattern Check (Used pattern)	Test pattern 1 Test pattern 2 Test pattern 3	Test pattern 1 Test pattern 2 Test pattern 3	Test pattern 1 Test pattern 2 Test pattern 3
Luminance Check	$L'max / L'min > 100$ $L'max > 200 \text{ cd/m}^2$ $\Delta(L'max / L'min) < 30 \%$ $\Delta Lamb < 30 \%$	$L'max / L'min > 40$ $L'max > 120 \text{ cd/m}^2$ $\Delta(L'max / L'min) < 30 \%$ $\Delta Lamb < 30 \%$	$L'max / L'min > 250$ $L'max > 250 \text{ cd/m}^2$ $L'min > 1.0 \text{ cd/m}^2$ $\Delta(L'max / L'min) < 30 \%$ $\Delta Lamb < 30 \%$
Grayscale Check	-	-	-
Uniformity Check	-	-	-
Multi-monitor	-	-	$\Delta L'max < 10 \%$ between multiple monitors $\Delta(L'max / L'min) < 10 \%$ between multiple monitors

Correlation between QS-RL and RadiCS

Pattern Check

The test patterns used are the same as the one specified in DIN V 6868-57.

Luminance Check

L_{max} and L_{min} used in QS-RL include the ambient luminance and are the same as L'_{max} and L'_{min} used in RadiCS.

QS-RL specifies $L_{min} \geq 1.0 \text{ cd/m}^2$, but RadiCS includes no equality sign. L_{amb} stands for the ambient luminance and refers to the same value as "L_s" of DIN V 6868-57. The inequality $L_{max}/L_s > 100$ (or 40) in the standard has been transformed into $L_s > L_{max}/100$ (or 40). In QS-RL, the luminance is specified as $|\Delta L_s| \leq 0.3 L_s$. Delta L_{amb} in RadiCS corresponds to the calculation of $|\Delta L_s|/L_s$ in QS-RL, and is expressed as its percentage. Note that none of inequalities used in RadiCS includes an equality sign.

L'_{max}/L'_{min} stands for a contrast ratio. The inequality of L'_{max}/L'_{min} in QS-RL has an equality sign in it ($L'_{max}/L'_{min} \geq 100, 40$ or 250) but the inequality in RadiCS does not. ($L'_{max}/L'_{min} > 100, 40$ or 250). In QS-RL, the luminance is specified as $|\Delta K_m| \leq 0.3 K_m$. K_m corresponds to L'_{max}/L'_{min} in RadiCS, and Delta (L'_{max}/L'_{min}) in RadiCS corresponds to the calculation of $|\Delta K_m|/K_m$ in QS-RL, and is expressed as its percentage. Also note that none of inequalities used in RadiCS includes an equality sign.

In QS-RL, L'_{max} and L'_{min} are determined by measuring the luminance at square regions filled with white (grayscale: 255) and black (grayscale: 0) in the test pattern 2, respectively. RadiCS displays 10 % of a display area in the middle and measures luminance by changing the grayscale 0 to 255. By doing so, the exact contrast ratio can be acquired.

Uniformity Check

The luminance uniformity is determined by firstly measuring the luminance of the center and a corner of the screen. Then, calculate the difference of these two luminance values and evaluate a percentage by dividing the difference by the luminance of the center. However, QS-RL does not specify particular measuring points for the uniformity measurement. In QS-RL, the measuring points are indicated with the test pattern 1 or the SMPTE pattern of the aspect ratio of 1:1, but the measuring points in these patterns have a significant difference, and other patterns around the measuring points may affect the measurement results. RadiCS displays two windows (grayscale: 128) with the size of 10 % of the whole display area at the center and a corner of the screen. It then measures the luminance at the center of both windows.

Since any monitors that support the RadiCS luminance check are LCD monitors, the criteria of 15 % or 20 % should apply to the LCD monitors, not to CRT monitors.

RadiCS specifies $(L_{corner}-L_{center})/L_{center} \times 100 < 15 \%$ (or 20 %), but this denotes $\pm 15 \%$ (or $\pm 20 \%$), and does not include an equals sign.

Sensors

DIN V 6868-57 requires the use of a measurement device for the acceptance tests that provides a luminance meter compliant with Class B or higher standard (DIN 5032-7) and does not block the ambient light. This requirement is also effective for QS-RL. RadiCS only allows noncontact type measurement devices to perform both the acceptance tests and the consistency tests. Since the EIZO sensors (UX2) are contact type measurement devices, they are not applicable.

Multi-monitor

Category A Mammo includes a determination for multiple monitors. QS-RL includes an equality sign but RadiCS does not.

Cautions

Category A Mammo conforms to the PAS1054 mammography standard. This standard includes monitor resolution of 2000 x 2500 or above as a condition, but RadiCS makes no such determination.

Basic QC

RadiCS Setup

	Acceptance Test	Consistency Test
Pattern Check (Used pattern)	TG18-QC	TG18-QC
Luminance Check	-	-
Grayscale Check	-	-
Uniformity Check	-	-
Multi-monitor	-	-

Correlation between Basic QC and RadiCS

Pattern Check

The patterns used for Basic QC are the same as those used for AAPM. RadiCS determines the properties to be verified and independently prepares appropriate patterns for each resolution.

- TG18-QC
This is scaled to match the resolution.

Basic Mammo QC and Basic Mammo QC for Remote

RadiCS Setup

	Acceptance Test	Consistency Test
Pattern Check (Used pattern) ^{*1}	TG18-QC TG18-UN80	TG18-QC TG18-UN80
Luminance Check	$L'max / L'min > 250$ $L'max > 450 \text{ cd/m}^2$ $Lamb < Lmin / 1.5$	$L'max / L'min > 250$ $L'max > 450 \text{ cd/m}^2$ $Lamb < Lmin / 1.5$
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 10 % of GSDF

^{*1} Not included in the consistency test of Basic Mammo QC for Remote.

Correlation between Basic Mammo QC, Basic Mammo QC for Remote and RadiCS

Pattern Check

The patterns used for Basic Mammo QC, Basic Mammo QC for Remote are the same as those used for ACR. RadiCS determines the properties to be verified and independently prepares appropriate patterns for each resolution.

- TG18-QC
This is scaled to match the resolution.
- TG18-UN80
A pattern solidly filled with white of grayscale 204.

Luminance Check

Except for $Lamb < Lmin / 1.5$, Basic Mammo QC and Basic Mammo QC for Remote include an equality sign in each judgment condition but RadiCS does not.

Grayscale Check

Basic Mammo QC and Basic Mammo QC for Remote include an equality sign but RadiCS does not because the target error rate is < 10 % of GSDF. This is a judgment condition for DICOM Part 14 GSDF. The number of grayscale measuring points is fixed at 18 and is unchangeable. The measurement result is 17 points because it is expressed as $(JND_{n+1} - JND_n)/2$.

Sensors

Any sensors can be used to perform both the acceptance tests and the consistency tests with Basic Mammo QC and Basic Mammo QC for Remote.

Basic QC Primary, Basic QC Primary for Remote, Basic QC Secondary, and Basic QC Secondary for Remote

RadiCS Setup

	Acceptance Test	
	Basic QC Primary Basic QC Primary for Remote	Basic QC Secondary Basic QC Secondary for Remote
Pattern Check (Used pattern)	-	-
Luminance Check	$L'_{\max} / L'_{\min} > 250$ $L'_{\max} > 170 \text{ cd/m}^2$ $\Lambda_{\text{mb}} < L_{\min} / 1.5$ $\Delta L'_{\max} < 10 \%$	$L'_{\max} / L'_{\min} > 100$ $L'_{\max} > 100 \text{ cd/m}^2$ $\Lambda_{\text{mb}} < L_{\min} / 1.5$ $\Delta L'_{\max} < 10 \%$
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 20 % of GSDF

	Visual Check ^{*1}	
	Basic QC Primary	Basic QC Secondary
Pattern Check (Used pattern)	TG18-QC	TG18-QC

	Consistency Test	
	Basic QC Primary Basic QC Primary for Remote	Basic QC Secondary Basic QC Secondary for Remote
Pattern Check (Used pattern) ^{*1}	TG18-QC	TG18-QC
Luminance Check	$L'_{\max} / L'_{\min} > 250$ $L'_{\max} > 170 \text{ cd/m}^2$ $\Lambda_{\text{mb}} < L_{\min} / 1.5$ $\Delta L'_{\max} < 10 \%$	$L'_{\max} / L'_{\min} > 100$ $L'_{\max} > 100 \text{ cd/m}^2$ $\Lambda_{\text{mb}} < L_{\min} / 1.5$ $\Delta L'_{\max} < 10 \%$
Grayscale Check	Target error rate < 10 % of GSDF	Target error rate < 20 % of GSDF

^{*1} Not included in Basic QC Primary for Remote and Basic QC Secondary for Remote.

Correlation between Basic QC Primary, Basic QC Primary for Remote, Basic QC Secondary, and Basic QC Secondary for Remote and RadiCS

Pattern Check

The patterns used for Basic QC Primary and Basic QC Secondary are the same as those used for AAPM.

RadiCS determines the properties to be verified and independently prepares appropriate patterns for each resolution.

- TG18-QC

This is scaled to match the resolution.

Luminance Check

Except for $\Lambda_{\text{mb}} < L_{\min} / 1.5$, Basic QC Primary, Basic QC Primary for Remote, Basic QC Secondary, and Basic QC Secondary for Remote include an equality sign in each judgment condition but RadiCS does not.

Grayscale Check

Basic QC Primary, Basic QC Primary for Remote, Basic QC Secondary, and Basic QC Secondary for Remote includes an equality sign but RadiCS does not because the target error rate is $< 10\%$ of GSDF. This is a judgment condition for DICOM Part 14 GSDF. The number of grayscale measuring points is fixed at 18 and is unchangeable. The measurement result is 17 points because it is expressed as $(JND_{n+1} - JND_n)/2$.

Sensors

Any sensors can be used to perform both the acceptance tests and the consistency tests with Basic QC Primary, Basic QC Primary for Remote, Basic QC Secondary, and Basic QC Secondary for Remote.

About Pathology350, Pathology450

RadiCS Setup

	Acceptance Test	
	Pathology350	Pathology450
Luminance Check	$L'_{\max} / L'_{\min} > 250$ $L'_{\max} > 350 \text{ cd/m}^2$ $\Delta L'_{\max} < 10 \%$	$L'_{\max} / L'_{\min} > 250$ $L'_{\max} > 450 \text{ cd/m}^2$ $\Delta L'_{\max} < 10 \%$
Grayscale Check	Target error rate < 10 %	Target error rate < 10 %

	Consistency Test	
	Pathology350	Pathology450
Luminance Check	$L'_{\max} / L'_{\min} > 250$ $L'_{\max} > 350 \text{ cd/m}^2$ $\Delta L'_{\max} < 10 \%$	$L'_{\max} / L'_{\min} > 250$ $L'_{\max} > 450 \text{ cd/m}^2$ $\Delta L'_{\max} < 10 \%$
Grayscale Check	Target error rate < 10 %	Target error rate < 10 %

Correlation between Pathology350, Pathology450 and RadiCS

Luminance Check

Pathology350, Pathology450 includes an equality sign in each judgment condition but RadiCS does not.

Grayscale Check

Pathology350, Pathology450 includes an equality sign but RadiCS doesn't because of the target error rate is < 10 %. The number of grayscale measuring points is fixed at 18 and is unchangeable.

Sensors

For Pathology350, Pathology450, any sensor can be used to perform both acceptance tests and consistency tests.

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Copyright © 2022 - 2025 EIZO Corporation. All rights reserved.
Java, Oracle Corporation ve/veya bağlı kuruluşlarının tescilli ticari markasıdır.
Intel, Intel Corporation şirketinin Amerika Birleşik Devletleri'nde ve/veya diğer ülkelerde ticari markasıdır.
CD mon, PEHA med Geräte GmbH şirketinin tescilli ticari markasıdır.
RaySafe, Unfors RaySafe AB şirketinin tescilli ticari markasıdır.
KONICA MINOLTA, Konica Minolta, Inc. şirketinin tescilli ticari markasıdır.