

## VisionCore FCS Viewer Supported Formats

### ● File Format for Loading (Video)

Video	Audio	Container	
H.264	MPEG Audio Layer 2, AAC	TS, MP4	
HEVC (H.265)	MPEG Audio Layer 2, AAC	TS, MP4	
Windows Media Video	Windows Media Audio	Windows Media	
Video for Windows	PCM	AVI	*1

\*1 Requires installation of plugins and codecs.

### ● File Format for Output (Video)

Video	Audio	Container	Mask, Annotation superimposed output	
H.264	MPEG Audio Layer 2, AAC	TS, MP4	Available	*1
HEVC (H.265)	MPEG Audio Layer 2, AAC	TS, MP4	Available	*2
MPEG-TS non-conversion cutout	-	TS	Unavailable	
Container conversion from MPEG-TS to MP4	-	MP4	Unavailable	
Windows Media Video	Windows Media Audio	Windows Media	-	
Video for Windows	PCM	AVI	-	*3

\*1 Supports Intel® Quick Sync Video (QSV), NVIDIA® NVENC, and software encoding

\*2 Supports Intel® Quick Sync Video (QSV) and NVIDIA® NVENC

\*3 Requires installation of plugins and codecs

### ● File Format for Loading (Still Image)

File extensions
.png .bmp .jpeg .tiff

### ● File Format for Output (Still Image)

File extensions
.png .jpeg .jpeg(4:2:2) .tiff

