



EIZO ColorEdge CG211 Wins TIPA 2007 Award for Best Monitor

Hakusan, Japan, May 2, 2007 – Eizo Nanao Corporation (TSE: 6737) is pleased to announce that its ColorEdge CG211 LCD monitor has been voted as the best monitor by TIPA (Technical Image Press Association) in its Photo and Imaging Awards 2007.

TIPA is a non-profit organization comprised of 31 European photo and imaging magazines with a readership of approximately eight million people in 12 countries. TIPA editors meet once a year to vote for the best photo and imaging products introduced in Europe in the past 12 months.

The ColorEdge CG211 is a 21.1-inch monitor for color work such as soft proofing and high-end digital photography. It features 12-bit hardware calibration capability and uniform brightness and color across the screen. EIZO's ColorNavigator calibration software is bundled with the CG211 and can be used with X-Rite (including GretagMacbeth and Monaco) and ColorVision measurement devices to calibrate the monitor in less than five minutes. The CG211 also offers hardware-based support for portrait mode with both Macintosh and Windows PCs.

About EIZO

Eizo Nanao Corporation is a leading global manufacturer of high-end visual display products with a wide range of LCD monitor systems. The image quality, long-term reliability, and innovative features of EIZO monitors make them the products of choice in many financial trading rooms, hospitals, back offices, and design studios. EIZO is based in Japan and represented in more than fifty countries by a network of exclusive distributors.

For more information, please contact:

Eizo Nanao Corporation
153 Shimokashiwano
Hakusan, Ishikawa 924-8566
Japan
Phone: +81 76 277-6792
Fax: +81 76 277-6793
<http://www.eizo.com/>

EIZO and ColorEdge are registered trademarks of Eizo Nanao Corporation.